AGU[®] ADVANCED
GROWING
UP

Balloon

NEBULIZER

COMPRESSOR

AGU N3

DE

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FR

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EN

USER MANUAL
Version 1.0

RU

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CONTENTS

| | | |
|----|--|----|
| 1 | Introduction..... | 45 |
| 2 | List of symbols | 45 |
| 3 | Field of application | 46 |
| 4 | Guidelines for safe operation | 46 |
| 5 | Complete set. Service life..... | 48 |
| 6 | Preparation for use | 50 |
| 7 | Operating principle and procedure..... | 51 |
| 8 | Adapter for mains..... | 52 |
| 9 | Cleaning..... | 53 |
| 10 | Filter change | 54 |
| 11 | Troubleshooting..... | 55 |
| 12 | Main technical characteristics..... | 56 |
| 13 | Technical details | 57 |
| 14 | Warranty..... | 62 |

1 INTRODUCTION

Thank you for purchasing Compressor nebulizer **AGU N3**. Proper use and care of your nebulizer will ensure effective and reliable medical treatment for many years.

2 LIST OF SYMBOLS

Symbol Meaning

The product complies with Directive 93/42/EEC on Medical Devices.

The product complies with TR CU 020/2011 «Electromagnetic compatibility of technical equipment».

WEEE (Waste Electrical and Electronic Equipment Directive). The symbol on the product or its package means that this product does not fall under the category of domestic waste. To avoid possible damage to the environment and human health, separate such wastes from others and dispose of them in accordance with accepted standards.

Type-BF device.

IP21

Ingress protection rating.

Leading digit (protection against ingress of solid foreign objects): 2 – protection against ingress of solid objects more than 12 mm in size; fingers or other objects with a maximum length of 80 mm, or solid objects.

Second digit (protection against vertically falling water drops): 1 – Vertically falling water drops should not impair the operation of the device.

Device class II (electric shock protection class).

Note/Warning.

Keep away from direct sunlight.

Read and understand the User Manual before use.

European Authorized Representative.

Manufacturer.

3 FIELD OF APPLICATION

Compressor nebulizer **AGU N3** is a medical device. Please use the device only under your doctor's and/or pulmonologist's direction. It provides comfortable, rapid and safe treatment. This device is suitable for all ages. Please read carefully this manual to familiarize yourself with features of the nebulizer.

NOTE

The nebulizer is designed for inhalation therapy of asthma, chronic obstructive pulmonary disease (COPD) and other respiratory tract diseases. Please consult your doctor before using the nebulizer for inhalation therapy. Follow the doctor's instructions regarding a dose and schedule of a medicine inhaled administration. This device corresponds to the Directive EC 93/42/EEC on Medical Devices, and BS EN 13544-1:2007+A 1:2009 Respiratory equipment – Part 1: Nebulizing systems and their components.

Please read carefully this manual before using the nebulizer.

NOTE

Please observe safety measures when using the nebulizer. The device must be used only for its intended purpose and with the prescribed medicines; according to this manual and in accordance with a doctor's instructions. This device is not intended for artificial lung ventilation and inhalation anaesthesia.

4 GUIDELINES FOR SAFE OPERATION

- Protect the device against water ingress to prevent electrical shock.
- Do not immerse the device in liquid.
- Do not use while bathing.
- Do not touch the device that has fallen into water. Unplug immediately.
- Do not use the previously immersed in water or dropped device with damaged parts (including power cable and plug). Please contact your service center to eliminate failures.

- Do not use the device when a malfunction is detected until the device is repaired.
- Do not use the device in areas with presence of flammable gases, oxygen or aerosol products.
- Do not use the device with empty medicine cup.
- Do not block the air openings. Do not place the device on a soft surface as the air openings may be blocked.
- Unplug the device before cleaning, filling and after each use.
- Unplug the device when not in use.
- Do not turn down or shake the device during operation.
- Use only accessories recommended by the manufacturer.
- Do not touch the power switch with wet hands.
- Do not disassemble or attempt to repair the device yourself.
- Do not use this device for artificial lung ventilation and inhalation anaesthesia.

USAGE PRECAUTIONS

- This device is not intended for use by people (including children) with physical, neurological or mental disorders, or who have lack of experience and knowledge, except when these persons are supervised or instructed on the use of the device by the person responsible for their safety. It is necessary to supervise children in order to prevent any improper use of with the device, its accessories, as well as factory package.
- Avoid contact of medicines, vapours or aerosol with eyes.
- The maximum capacity of the medicine cup is 6 ml. Do not exceed this amount.
- Do not use the device while driving.
- Unplug the device as soon as you feel discomfort or experience reduced general condition during inhalation.
- Do not use the device if the air hose is curved.

STORAGE PRECAUTIONS

- Keep away from direct sunlight, high temperatures and humidity.
- Keep out of the reach of children.
- Unplug the device when not in use.

CLEANING PRECAUTIONS

- Check air filter, mouthpiece, nebulizer and other auxiliary components before each use. Contaminated or worn out components must be replaced.
- Do not immerse the device in liquid to avoid its damage.
- Unplug the device before cleaning.
- Clean all necessary accessories after each use according to the instructions in this manual.

- Remove any medicine residues from the medicine cup after each use. Always check the expiration date of the medicines before inhalation.
- Do not store wet air hose or air hose with medicine residues. This may cause bacterial infection.

MEDICAL PRECAUTION

This manual should not replace a consultation by a medical specialist. The data contained in this document, or this device should not be used for self-diagnosis and treatment, or as a basis for choice of a medicine. At suspicion on any disease requiring inhalation of a medicine, please consult a doctor.

5 COMPLETE SET. SERVICE LIFE

The complete set includes the following components. If you discovered any missing components, please contact the seller immediately.

The nebulizer system consists of the following components:

- Compressor (main unit) – 1 pc.
- Mains adapter – 1 pc.
- Power cable – 1 pc.

Periodically replaced parts and components:*


- Children mask – 1 pc.
- Adult mask – 1 pc.
- Nebulizer – 1 pc.
- Baffle – 1 pc.
- Nebulizer cover – 1 pc.
- Mouthpiece – 1 pc.
- Nosepiece – 1 pc.
- Air hose – 1 pc.
- Replaceable air filters – 5 pcs.
- Mains adapter – 1 pc.
- Nozzle – 1 pc.

*Parts and components replacement interval is determined by their service life specified in the section below.

The service life of the components may vary depending on the intensity of use. When using the device for spraying 2 ml of saline 2 times a day for 10 minutes at room temperature (23 °C), the service life of the components is as follows:

- | | |
|--|----------|
| • main unit | 3 years |
| • nebulizer, baffle, nebulizer cover, nozzle | 1 year |
| • air hose, mouthpiece, nosepiece | 1 year |
| • air filter | 30 days* |
| • adult/children mask | 1 year |
| • mains adapter | 1 year |

*Note: the air filter service life may be significantly reduced when using the device in an environment with increased dust concentration.

6


PREPARATION FOR USE

Before first use or after long-term non-use of the nebulizer it must be prepared for use according to the instructions set in the «Cleaning» section of this manual.

NOTE

Before cleaning, assembling and after each use of the device make sure it is turned off and unplugged.



- 1 Place the compressor on a flat, stable surface.
- 2 Gently open the nebulizer cover to separate it into two sections (medicine cup and cover).
- 3 Be sure the baffle is properly inserted into the medicine cup. The stem inside the medicine cup must be inserted into the groove of the baffle.
- 4 Gently twist the medicine cup and the cover together to assemble the nebulizer. Be sure the two sections are fit securely.

7

OPERATION PRINCIPLE AND PROCEDURE

The Compressor nebulizer may be operated at an inclination angle of up to 45° against the vertical line. If the inclination angle exceeds 45°, aerosol will not be generated.

- 1 Add the inhalation solution to the medicine cup. Capacity of the medicine cup is 6 ml.
- 2 Connect one end of the air hose to the air hose connector.
- 3 Connect the other end of the air hose to the compressor opening.
- 4 Attach the mouthpiece, mask or nosepiece to the upper part of the nebulizer.


NOTE

If several inhalation sessions are required, it is recommended to make a 30-minute break after each session.

After each use:

- turn off the device and unplug the power cable;
- let the device cool down completely;
- gently disconnect the air hose from the nebulizer and empty the residues of the medicine;
- clean the device according to the instructions of this manual.

8 MAINS ADAPTER


9

CLEANING

FLUSHING

Flush the device after each procedure or before its first use.

- 1 Make sure the device is turned off and unplugged.
- 2 Disconnect the air hose from the base of the nebulizer.
- 3 Gently twist and pull the nebulizer cover to open it.
- 4 Rinse the nebulizer sections under running hot water.
- 5 Wipe the sections with a clean towel or allow them to dry completely in the air.
- 6 Assemble the nebulizer.

NOTE

Clean thoroughly all the nebulizer sections and pieces/caps, including the air hose, during the first cleaning and after a long-term storage of the device.

CLEANING OF COMPRESSOR


Wipe the dirty compressor with a soft cloth.

NOTE

Any other cleaning method or use of a detergent may affect the appearance and performance of the device.

10 FILTER CHANGE

The air filter must be replaced when it turns grey.



1 Gently pull the air filter cover to remove it.

2 Remove the dirty filter.

3 Insert new, clean filter.

4 Close tight the air filter cover.

NOTE

Do not wash or clean the air filters. Only **AGU N3** air filters can be used. Do not use alternative filter materials like cotton cloth. Do not operate the device without the air filter.

11 TROUBLESHOOTING

Please note that some problems can be solved without a specialist involvement.

| Problem | Solution |
|---|---|
| The device does not operate as it is not plugged into the mains | Make sure the adapted is plugged into the socket |
| The device does not nebulize or it nebulizes weakly | Check the medicine in the cup |
| | Check the nebulizer for damage |
| | Check the position of the baffle inside the nebulizer |
| | Make sure the air hose and other auxiliary accessories are properly connected |
| | Check the air filter and replace if necessary |

Electric shock protection class:

Device class II.

Type-BF device:

Mouthpiece and mask.

Ingress protection rating:

IP21.

Degree of safety in the presence of flammable aesthetic agents or oxygen:

Not AP/APG (do not use with flammable aesthetic agents or oxygen).

12 MAIN TECHNICAL CHARACTERISTICS

| Type | Compressor nebulizer |
|------------------------------------|---|
| Model | AGU N3 |
| Mains adapter | DC 5 V, 2 A |
| Voltage | AC 100~240 V, 50-60 Hz |
| Wattage, W | ≤ 5 |
| Sound level (1 m from AGU N3), dBA | ≤ 45 |
| Air flow rate, l/min | ≥ 2.5 |
| Operating temperature | +10 °C – +40 °C |
| Humidity | ≤ 90% RH |
| Operating atmosphere pressure, hPa | 700–1060 |
| Storage conditions | –10 °C – +60 °C, permitted RH ≤ 90% |
| Dimension (LxWxH), mm | 98x57x50 |
| Weight (without accessories), g | 170 |
| Capacity of medicine cup, ml | 6 |
| Particle size (MMAD), μm | ≤ 2.9 |
| Average nebulizing rate, ml/min | ≥ 0.3 (0.9 % of NaCl) |
| Complete set | Nebulizer, air hose, mouthpiece, nosepiece, air filters (5 pcs.), children mask, adult mask, mains adapter, nozzle, bag |
| Service life | Main unit – 3 years; nebulizer, baffle, nebulizer cover, nozzle – 1 year; air hose, mouthpiece, nosepiece – 1 year; air filter – 30 days*; adult/children mask – 1 year; mains adapter – 1 year |

***Note:** the air filter service life may be significantly reduced when using the device in an environment with increased dust concentration.

The manufacturer reserves the right for technical changes in the design and the complete set of the device.


Depending on a medicine (aqueous solution, suspension or high-viscosity substance), the aerosol delivery may vary.

13 TECHNICAL DETAILS

Measurements of the particle size using the cascade impactor for **AGU N3** according to BS EN 13544-2009 «Respiratory equipment – Part 1: Nebulizing systems and their components», Annex CC.3, based on the measurements of the particle size with the cascade impactor are as follows:

| | |
|------------------------------------|--|
| Size of nebulized particles | MMAD $\leq 2.9 \mu\text{m}$ (MMAD – Mass Median Aerodynamic Diameter) |
| Capacity of medicine cup | Max. 6 ml |

Distribution of particles

ELECTROMAGNETIC EMISSIONS: GUIDANCE AND MANUFACTURER'S DECLARATION

Compressor nebulizer **AGU N3** is intended for use in the electromagnetic environment specified below.

The customer or user of **AGU N3** should assure that it is used in such environment.

| Emission test | Compliance | Electromagnetic environment – Guidance |
|---|------------|---|
| RF emissions CISPR 11 | Group 1 | Compressor nebulizer AGU N3 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment |
| RF emissions CISPR 11 | Class B | Compressor nebulizer AGU N3 is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes |
| Harmonic emissions IEC 61000-3-2 | Class A | |
| Voltage fluctuations/ flicker emissions IEC 61000-3-3 | Complies | |

| Immunity test | IEC 60601 Test level | Compliance level | Electromagnetic environment – guidance |
|--|---|---|---|
| Electrostatic discharge (ESD) IEC 61000-4-2 | ± 6 kV contact ± 8 kV air | ± 6 kV contact ± 8 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 % |
| Electrical fast transient/burst IEC 61000-4-4 | ± 2 kV for power supply lines ± 1 kV for input/output lines | ± 2 kV for power supply lines Not applied | Mains power quality should correspond to the standards of the public power supply networks |
| Surge IEC 61000-4-5 | ± 1 kV line(s) to line(s) ± 2 kV line(s) to earth | ± 1 kV line(s) to line(s) Not applied | Mains power quality should correspond to the standards of the public power supply networks |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | <5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 s | <5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 s | Mains power quality should correspond to the standards of the public power supply networks. If the user of the device requires continued operation during power mains interruptions, it is recommended that AGU N3 to be powered from an uninterruptible power supply or a battery |
| Power frequency magnetic fields (50/60 Hz) IEC 61000-4-8 | 3 A/m | 3 A/m | Power frequency magnetic fields of AGU N3 should be at levels characteristic of a typical public use systems |

NOTE: U_T – is the AC mains voltage prior to application of the test level.

| Immunity test | IEC 60601 Test level | Compliance level | Electromagnetic environment – guidance |
|---|---|---------------------|--|
| Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3 | 3 Vrms 150 kHz – 80 MHz 3 V/m 80 MHz – 2.5 GHz | 3 Vrms 3 V/m | <p>Portable and mobile RF communications equipment should be used no closer to any part of AGU N3 including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance: $d = 1.2 \sqrt{P}$, $d = 1.2 \sqrt{P}$ 80 MHz – 800 MHz, $d = 2.3 \sqrt{P}$ 800 MHz – 2.5 GHz, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters as determined by an electromagnetic site survey a should be less than the compliance level in each frequency range b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p> |

NOTE 1: at 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: these guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

- a.** Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which **AGU N3** is used exceeds the applicable RF compliance level above, **AGU N3** should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating **AGU N3**.
- b.** Over the frequency range 150 kHz to 80MHz, field strengths should be less than 3 V/m.

RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND COMPRESSOR NEBULIZER AGU N3

Compressor nebulizer **AGU N3** is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of compressor nebulizer **AGU N3** can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and compressor nebulizer **AGU N3** as recommended below, according to the maximum output power of the communications equipment. For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

| Rated maximum output power of transmitter (W) | Separation distance according to frequency of transmitter (m) | | |
|---|---|--|---|
| | 150 kHz – 80 MHz $d = 1.2 \sqrt{P}$ | 80 MHz – 800 MHz $d = 1.2 \sqrt{P}$ | 800 MHz – 2.5 GHz $d = 2.3 \sqrt{P}$ |
| 0.01 | 0.12 | 0.12 | 0.23 |
| 0.1 | 0.38 | 0.38 | 0.73 |
| 1 | 1.2 | 1.2 | 2.3 |
| 10 | 3.8 | 3.8 | 7.3 |
| 100 | 12 | 12 | 23 |

NOTE 1: at 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: these guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

WARNING

This symbol on the product means that this product is an electronic device which must be disposed of in accordance with Directive 2012/19/EU at the local waste disposal utility.

14 WARRANTY

The warranty period for the main unit of compressor nebulizer **AGU N3** is 24 months from the date of sale, providing that all operation conditions specified in this instruction manual are strictly observed. The warranty is only valid if there is a cashier's receipt and a warranty certificate completed by an official representative to confirm the date of sale. The warranty does not apply to the consumables, parts and components subject to periodic replacement, such as: children mask, adult mask, replaceable filters, nebulizer, air hose, mouthpiece, nosepiece, package.

After-sales and free maintenance service is not provided in case of:

- using the device with violation of requirements specified in the instruction manual;
- damage caused by deliberate or erroneous actions of the consumer due to mishandling or negligence;
- evidence of mechanical damage, dents, cracks, chips, etc. on the device casing, opening of the casing, disassembly, unauthorized repair, ingress of moisture into the casing or effect of corrosive substances, or any other unauthorized interference, as well as in other cases of violation by the consumer of storage, cleaning, transportation and operation requirements specified in the instruction manual;
- ingress of oils, dust, insects, liquids (not intended for use with the device) and other foreign objects into the device.

The warranty does not cover the defects or malfunction caused by:

- normal wear and tear of the components with limited service life;
- damage of accessories which are not an integral part of the device (children mask, adult mask, replaceable filters, nebulizer, air hose, mouthpiece, nosepiece, package);
- using defective, worn out accessories and pieces;
- force majeure circumstances (accident, fire, flood, electric line fault, etc.).

Date of manufacture  and importer information are specified on the individual package.

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AGU[®] is the registered trademark by Montex Swiss AG,
Donawa Lifescience Consulting Srl,