Brainy

SMART THERMOMETER

INFRARED

AGU SHE10

de

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CONTENTS

1	Introduction.....	32
2	List of symbols	33
3	Scope of application	35
4	Complete set	35
5	Basic specifications	36
6	Device description	37
7	Guidelines for safe operation	39
8	Before you start	40
9	Operation principle and procedure	41
10	Cleaning	49
11	Troubleshooting	50
12	Storage, operation and transportation rules	51
13	Disposal	52
14	Certification	52
15	Manufacturer's warranty	56

1 INTRODUCTION

Dear friends, thank you for choosing our products!

Smart infrared thermometer **AGU SHE10** has a wide measurement range, which allows you to use it both for measuring body temperature (in the ear or on the forehead) and for measuring:

- milk surface temperature in a baby bottle;
- water surface temperature in a baby bathtub;
- ambient temperature.

Innovative infrared sensor technology allows for measuring temperature on the forehead in a few seconds.

The unique design of the device with a built-in infrared sensor provides accurate and reliable measurement results.

It is easy and simple to measure temperature with Smart infrared thermometer **AGU SHE10**.

This thermometer has the ability to synchronize with a mobile application. This allows you to store measurement results and add notes to measurements.

Please read this user manual carefully before using the device.

2 LIST OF SYMBOLS

Symbol	Meaning
	CE marking with Notified Body identification number. Indicates compliance with the European Medical Device Directive 93/42/EEC.
	The product complies with the basic requirements of TR CU 020/2011 «Electromagnetic compatibility of technical means».
	WEEE (Waste Electrical and Electronic Equipment Directive). The symbol on the product or its package means that this product does not fall under the category of domestic waste. To avoid possible damage to the environment and human health, separate such wastes from others and dispose of them in accordance with accepted standards.
	Special disposal of batteries. The symbol on the product or its package indicates that the batteries are not classified as household waste. To avoid possible damage to the environment and human health, separate such wastes from others and dispose of them in accordance with accepted standards.
	Type BF device (body floating) – Type B device with an insulated F-Type applied part. Type of protection of applied part against electric shock.

IP22 Ingress protection rating. Leading digit (protection against ingress of solid foreign objects): 2 – protection against ingress of solid objects more than 12 mm in size; fingers or other objects with a maximum length of 80 mm, or solid objects. Second digit (protection against ingress of foreign liquids): 2 – protection against vertically dripping water at an angle of up to 15° from the vertical (equipment in normal position).

Warning/Attention.

Please read this user manual carefully before using the device.

Manufacturer.

Authorized manufacturer's representative in the European Union.

Temperature limits.

3 SCOPE OF APPLICATION


Smart infrared thermometer **AGU SHE10** is intended for domestic and clinical use. With the device you can measure body temperature both on the forehead and in the external ear canal. This thermometer can be used both for children and adults. Do not use the device for purposes other than intended.

4 COMPLETE SET

- 1 Smart infrared thermometer **AGU SHE10** – 1 pc.
- 2 1.5 V AAA battery – 2 pcs.
- 3 User manual with quick guide – 1 pc.

2 X 1.5 V AAA

5 BASIC SPECIFICATIONS

Type	Infrared digital thermometer
Model	AGU SHE10
Power supply	3 V (2 x 1.5 V AAA batteries)
Measurement range	Body temperature when measured on the forehead: 34.0 °C ~ 43 °C (93.2 °F ~ 109.4 °F) Object temperature: 2 °C – 93.2 °C (35.6 °F ~ 199.9 °F) Body temperature when measured in the external ear canal: 34.0 °C ~ 43 °C (93.2 °F ~ 109.4 °F)
Measuring accuracy	Objects: $T \leq 40 \text{ °C} \pm 2 \text{ °C}$ ($\pm 0.4 \text{ °F}$); $T > 40 \text{ °C} \pm 5$ Human body when measured on the forehead: $\pm 0.2 \text{ °C}$ ($\pm 0.4 \text{ °F}$) Human body when measured in the external ear canal: $\pm 0.2 \text{ °C}$ ($\pm 0.4 \text{ °F}$)
Dimensions, mm	~145 (L) x ~38 (W) x ~50 (H)
Weight (with batteries), g	Approx. 100
Housing material	ABS plastic, acryl
Measurement time	From 1.4 sec
Display	Segment display with special icons
Auto power-off	30 seconds
Operating conditions	Temperature: 15 °C ~35 °C (59 °F ~ 95 °F) Humidity: up to 95% (non-condensing)
Memory	9 measurements

Indication when measuring body temperature

Temperature	Indication
$T < 34.0 \text{ °C}$ (93.2 °F)	«L»
34.0 °C (93.2 °F) $\leq T \leq 43.0 \text{ °C}$ (109.4 °F)	Shows temperature value
$T > 43.0 \text{ °C}$ (109.4 °F)	«H»

Acoustic signals


The thermometer is switched on and ready for use	1 short beep
Measurement completed	5 beeps if the value is above 37.5 °C
	1 long beep if the value is below 37.5 °C
System error or malfunction	3 short beeps

6 DEVICE DESCRIPTION

Smart infrared thermometer **AGU SHE10** features:


- temperature measurement within 1 second;
- high measurement accuracy;
- sound indication.

Wireless synchronization with AGU app.

- 1 Cover cap.
- 2 IR sensor.
- 3 Battery compartment cover.
- 4 Segment display.
- 5 **START** button.
- 6 **ON/OFF** button.
- 7 Measurement mode switch.

DESCRIPTION OF DISPLAY ELEMENTS

- 1 Measurement result in Fahrenheit.
- 2 Measurement result in Celsius.
- 3 Low battery level.
- 4 External ear canal measurement mode.
- 5 Forehead measurement mode.
- 6 Measurement in process.
- 7 Object temperature measurement mode.
- 8 Sound notification enabled.
- 9 Measurement result.

7 GUIDELINES FOR SAFE OPERATION

Follow the instructions precisely to ensure reliable and long-term operation of the device.

- Any modification of the device is prohibited.
- Do not expose the thermometer to mechanical shock or use if damaged.
- Do not immerse the device in water or other liquids.
- Do not use the device near strong electromagnetic sources, such as wireless communication sources or cellular phones.
- Do not allow children to play with the thermometer.
- Do not expose the thermometer to direct sunlight, extreme temperatures, moisture or contamination.
- The environmental temperature in the apartment may differ in various rooms. Before measuring temperature, the user and the thermometer should stay in the same room (premises) with constant environmental conditions for at least 15 minutes.
- Before starting measurement, please remove any dirt, hair or sweat from the measuring area.
- Do not measure temperature during or immediately after breastfeeding.
- The user should not take food or drink immediately before or during temperature measurement.
- Do not pocket the thermometer before use as it can be heated.
- Every time before starting measurement, please make sure that the infrared sensor lens is not contaminated. Clean the measuring sensor lens if it is contaminated.


WARNING

During early stages of illness with a rapid rise in temperature to high values, the effect of «white hyperthermia» can be observed – a condition which can be characterized by the peripheral vessel constriction, when the skin becomes pale and remains cold. In such cases, temperature measurements on the forehead should not be carried out, because the skin temperature will be low.

8 BEFORE YOU START

- 1 Remove the device from the package.
- 2 Remove the battery compartment cover.
- 3 Insert new 2 x 1.5 V AAA batteries into the battery compartment, observing polarity. Always use only new batteries.
- 4 Close the battery compartment cover.

Dispose of used batteries in accordance with applicable legal regulations. It is forbidden to dispose of batteries in household waste.

9 DEVICE SETTINGS

DATE AND TIME SETTINGS

Before the first use of the thermometer and after each replacement of batteries, date and time are displayed in the format 1-1 00:00.

Perform a step-by-step setup of the device:

- Press and hold the **ON\OFF** button for 5 seconds.
- Set the year (from 20 to 35) by pressing the **START** button.
Press the **ON\OFF** button again.
- Set the month (from 1 to 12) by pressing the **START** button.
Press the **ON\OFF** button again.
- Set the day (from 1 to 31) by pressing the **START** button.
Press the **ON\OFF** button again.
- Set the hour (from 0 to 23) by pressing the **START** button.
Press the **ON\OFF** button again.
- Set the minutes (from 0 to 59) by pressing the **START** button.
Press the **ON\OFF** button again to complete the setup.

REPRODUCTION OF THE PREVIOUS READINGS

You can recall up to 9 measurement values.

- 1 Press the **START** button.
- 2 The display will show the memory location number, then the value.
- 3 Press the **START** button to switch memory locations.

When recording new measurement values (after the 9th value), old values are deleted automatically.

SOUND ON/OFF

- 1 Press and hold the **START** button for 5 seconds.
- 2 The display will show the beeper icon and its status.


Press the **START** button to turn the sound notification **ON** or **OFF**.

- 3 Do not press any buttons on the device for 3 seconds. The thermometer will turn off after saving the settings.

SWITCHING BETWEEN THE UNITS OF MEASUREMENT

- 1 Press and hold the **START** button for 9 seconds.
- 2 The display will show °C or °F.

- 3 Press the **START** button to switch between the units of measurement.
- 4 Do not press any buttons on the device for 3 seconds. The thermometer will turn off after saving the settings.

OPERATION PRINCIPLE AND PROCEDURE

There are 2 options for the device operation:

- Option 1:** without mobile application.
- Option 2:** via mobile application.

DEVICE OPERATION WITHOUT MOBILE APPLICATION

Smart infrared thermometer **AGU SHE10** measures the body temperature on the forehead as well as the temperature of objects.


MEASURING BODY TEMPERATURE ON THE FOREHEAD

- 1 Turn on the thermometer by pressing the **ON/OFF** button.
- 2 Remove any dirt, hair or sweat before placing the temperature sensor in the measurement area.
- 3 Flip the toggle switch on the back of the thermometer to the forehead measurement mode   .
- 4 The display will show the forehead temperature measurement icon and the object temperature measurement icon.

The thermometer is ready for measurement.

- 5 Put the sensor tip with the cap to the temple and while holding the **START** button smoothly slide over the forehead in the direction of the other temple.

The arrow on the display will flash during measurement. Release the **START** button.

- 6 When the measurement is complete, you will hear a beep and the display will show the measured body temperature.


MEASURING OBJECT TEMPERATURE

- 1 Turn on the thermometer by pressing the **ON/OFF** button.
- 2 Flip the toggle switch on the back of the thermometer to the forehead measurement mode   .
- 3 The display will show the forehead temperature measurement icon and the object temperature measurement icon.

The thermometer is ready for measurement.

- 4 Put the sensor tip with the cap to the object and press the **START** button twice.

- 5 You will hear a beep and the display will show the measured temperature value.

TEMPERATURE MEASUREMENT IN THE EXTERNAL EAR CANAL

WARNING

- Some people produce different readings in their left and right ear. In order to avoid changes in readings, always measure a person's temperature in the same ear.
- The measurement must not be taken in an ear affected by inflammatory diseases, after possible ear injuries or in the healing period. In all of these cases, please consult your doctor.
- If you have been lying on one ear for some time, the temperature is slightly raised. Wait a little while or measure in the other ear.
- Accumulation of ear wax on the sensor can affect the measuring accuracy and result in the spread of infection among the thermometer users. It is extremely important for each measurement to observe the rules of hygiene and keep the sensor clean. To clean the sensor, follow the instructions in the «Device cleaning» section. Every time before starting measurement, please make sure that the measuring sensor lens is not contaminated.
- After cleaning the measuring sensor lens with a disinfectant solution, wait a few minutes before taking the next measurement for the sensor to return to its initial operating temperature.


The thermometer measures the infrared heat generated by ear drum and surrounding tissues. The energy is captured by the lens and converted into a temperature value.


Temperature measurement should be carried out by placing the tip in the external ear canal above the ear drum, and not over the surface of nearby

anatomical areas. Failure to do so may result in excessive readings.

Body temperature measurement in the external ear canal ensures accurate readings due to the same temperature of blood flowing towards the ear drum and the brain's temperature regulation center – hypothalamus. Therefore, changes in body temperature are faster and more accurately recorded in the external ear canal than in any other body area.

The thermometer measures the infrared heat generated by ear drum and surrounding tissues.

- 1 Turn on the thermometer by pressing the **ON/OFF** button.
- 2 Flip the toggle switch on the back of the thermometer to the external ear canal measurement mode   .
- 3 The display will show the external ear canal temperature measurement icon. 
- 4 Remove the cap.

The thermometer is ready for measurement.

- 5 Make sure the ear canal is clean. As the ear canal is slightly curved, you have to pull the ear slightly up and backwards before inserting the sensor tip.


THE CORRECT POSITION OF THE THERMOMETER SENSOR IN THE EAR CANAL

- **Children under 1 year of age:**

Place your child in a horizontal position and turn the head on one side so that the ear is directed upward. Slightly pull the ear back. Carefully insert the sensor tip and press the **START** button.



- **Children over 1 year of age:**

Slightly pull the ear up and back to align the ear canal. Carefully insert the sensor tip and press the **START** button.

The arrow on the display will flash during measurement.

- 6 You will hear a beep and the display will show the measured temperature value.

DEVICE OPERATION VIA MOBILE APPLICATION

AGU mobile application provides opportunity to:

- save the temperature measurement history for different users;
- take notes on symptoms and general state;
- take notes on medicines taken.

ATTENTION

The object temperature measurement history is not stored in the application.

Application installation: Android 6.0 or higher, iOS 13.0 or higher.


Option a: open the PlayMarket or AppStore, enter «AGU» in the search box and download the application.

Option b: for direct link to download the application, scan this QR code.

- 1 Download and run the application.
- 2 Register or log in through social media site.
- 3 Select device to connect.
- 4 Enter the required user data.
- 5 Switch on the device and press «Connect» button.

- 6 Follow the application instructions.


10 DEVICE CLEANING

WARNING

Do not use chemically active detergents to clean the device.





- 1 Clean the device with a soft, dry cloth. If necessary, the cloth can be slightly dampened.
- 2 If any dirt gets on the infrared sensor, the readings may differ from the actual ones, as this affects the measurement (usually the readings are too low). To eliminate this problem, it is enough to wipe the infrared sensor with a soft lint-free cloth (or cotton swab) moistened with disinfectant solution. After the cleaning procedure, the infrared sensor lens surface should be crystal clear. As the lens surface gets cool when the disinfectant solution evaporates, the thermometer should be left for several minutes at room temperature before use after cleaning.

- 3 Leave the device to air dry at room temperature.
- 4 Store the thermometer in a dry, clean place in accordance with the operating conditions.

11 TROUBLESHOOTING

If a problem occurs during operation, first of all, refer to the list of possible failures given below.

Indication	Cause	Solution
	<p>The measured temperature is higher: in body temperature measurement mode – $T > 43.0\text{ °C}$ (109.4 °F)</p>	<p>Use the thermometer only for the specified temperature ranges: If necessary, clean the sensor tip. In the event of a repeated error message, contact your retailer or Customer Service</p>
	<p>The measured temperature is higher: in object temperature measurement mode – $T > 93.2\text{ °C}$ $(\text{or } 199.9\text{ °F})$</p>	
	<p>The measured temperature is lower: in body temperature measurement mode – $T < 34.0\text{ °C}$ (93.2 °F)</p>	
	<p>The measured temperature is lower: in object temperature measurement mode – $T < 2.0\text{ °C}$ (35.6 °F)</p>	

	<p>1. Operating temperature is out of range between 15 °C – 35 °C (59 °F – 95 °F)</p> <p>2. Too frequent temperature measurements</p>	<p>1. Use the thermometer only for the specified temperature ranges.</p> <p>2. It is necessary to maintain a pause of 5 seconds between measurements.</p>
	<p>Very low battery</p>	<p>Replace batteries</p>

12 STORAGE, TRANSPORTATION AND OPERATION RULES

- The device should be stored at temperatures between -25 °C and $+55\text{ °C}$ (-13 °F ... $+131\text{ °F}$) with a maximum relative humidity of 95% (non-condensing).
- The device should be operated at temperatures between $+15\text{ °C}$ and $+35\text{ °C}$ (59 °F – 95 °F) with a relative humidity of up to 95% (non-condensing).
- The device should be transported at temperatures between -25 °C and $+55\text{ °C}$ (-13 °F ... $+131\text{ °F}$) with a maximum relative humidity of 85% (non-condensing).
- Do not expose the device to thermal shock.

ATTENTION

After transportation or storage at low temperatures, it is necessary to keep the device at room temperature for at least 2 hours before switching on.

13 DISPOSAL

The symbol on the product or its package means that this product does not fall under the category of domestic waste.

- With proper disposal of the product, you will help to prevent possible negative effects of the device on the environment and human health.
- For the protection of the environment, the device should not be disposed of with domestic (household) waste. Disposal must be carried out in accordance with local regulations.
- Device must be disposed of in accordance with the EU Directive 2012/19/EU WEEE (Waste Electrical and Electronic Equipment).

If you have any questions, please contact the local public utility responsible for waste disposal.

Dispose of used, fully discharged batteries in special containers, take to special waste collection centers or electronics stores. Users are legally obliged to ensure the disposal of batteries.

14 CERTIFICATION

This product complies with the provisions of Medical Devices Directive (MDD) 93/42/EEC.

EMC Tables

AGU SHE10 is intended for use in the electromagnetic environment specified below. The customer or the user of **AGU SHE10** must make sure that it is used in such an environment

Guidance and manufacturer's declaration –
Electromagnetic emissions

Phenomenon	Professional healthcare facility environment a)	HOME HEALTHCARE ENVIRONMENT a)
Conducted and radiated RF EMISSIONS	a)	CISPR 11 Group 1 Class B
Harmonic distortion	Not applicable	
Voltage fluctuations and flickering	Not applicable	
<p>a) The equipment is suitable for use in Home Health Environments and Professional Health Care Environments limited to patient rooms and respiratory treatment facilities in hospital or clinics. The more restrictive acceptance limits of Group 1 Class B (CISPR 11) have been considered and applied. The equipment is suitable for use in the mentioned environments when directly connected to the Public Mains Network.</p> <p>b) The test is not applicable in this environment unless the ME EQUIPMENT and ME SYSTEM used will be connected to the PUBLIC MAINS NETWORK and the power input is otherwise within the scope of the Basic EMC standard</p>		
Guidance and manufacturer's declaration – Electromagnetic immunity – Enclosure port		

Phenomenon	Basic EMC standard or test method	Immunity test levels	
		Professional health-care facility environment	HOME HEALTHCARE ENVIRONMENT
ELECTROSTATIC DISCHARGE	IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	
Radiated RF EM fields	IEC 61000-4-3	a)	10 V/m b) 80MH -27 GHz 80% AM at 1kHz

Proximity fields from RF wireless communications equipment	IEC 61000-4-3	<p>COMPLIANT</p> <p>NOTE: Further information about distances to be maintained between portable and mobile RF communications equipment (transmitters) and the AGU SHE10 can be requested from supplier using the contact information provided in this manual. However, it is advisable to keep the electro-mechanical aerosol equipment at an adequate distance of, at least, 0.5 m from mobile phones or other RF communications transmitters to minimize possible interference</p>
RATED power frequency magnetic fields	IEC 61000-4-8	30 A/m c) 50 Hz or 60 Hz
<p>a) The equipment is suitable for use in Home Health Environments and Professional Health Care Environments limited to patient rooms and respiratory treatment facilities in hospital or clinics. The more restrictive IMMUNITY acceptance limits have been considered and applied.</p> <p>b) Before modulation is applied.</p> <p>c) This test level assumes a minimum distance of at least 15 cm between the ME EQUIPMENT or ME SYSTEM and sources of power frequency magnetic fields</p>		

15 MANUFACTURER'S WARRANTY

Warranty information is specified in the warranty card. The warranty period starts from the date of sale, providing that all operation conditions specified in this User Manual are strictly observed. The warranty is valid only on presentation of the sales receipt and the warranty card completed by the authorized dealer. The warranty does not apply to consumables, constituent parts and components, such as battery, packaging, etc.

Free warranty service is not provided in case of.

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- use of the device in violation of the operating instructions;
- damage caused by intentional or erroneous actions of the consumer due to improper or negligent handling; evidence of mechanical impact, dents, cracks, chips, etc. on the thermometer housing, evidence of opening the device, disassembly, unauthorized repair, ingress of moisture, exposure to aggressive agents or any other unauthorized alterations of the device, and in other cases of violation of storage, cleaning, transportation and operation rules specified in the User Manual;
- ingress of oils, dust, insects, liquids and other foreign objects inside the device.

Follow the instructions precisely to ensure reliable and long-term operation of the device.

In case of abnormal operation of the device, please contact your retailer or Customer Service.

The manufacturer reserves the right to make alterations in the design of the device.

WARNING

Do not attempt to repair the device yourself – this will void the warranty.

The thermometer alterations without the manufacturer's permission will void the warranty.

For repair and maintenance, please contact a specialized after-sales service (see agu-baby.com).

Date of manufacture and importer information are specified on the individual package.

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AGU[®] is the registered trademark by Montex Swiss AG,
AViTA Corporation, 9F, No. 78, Sec. 1, Kwang-Fu Rd.,