Giraffe

THERMOMETER

NON-CONTACT
FOR CHILDREN

AGU NC8/NC8C/NC8A

DE

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agu-baby.com

CONTENTS

1	Introduction.....	25
2	List of symbols.....	26
3	Scope of application	27
4	Complete set.....	27
5	Basic specifications	28
6	Device description.....	29
7	Guidelines for safe operation	32
8	Before you start.....	33
9	Operation principle and procedure	33
10	Device cleaning.....	39
11	Troubleshooting.....	40
12	Storage, transportation and operation rules.....	41
13	Disposal.....	41
14	Certification	42
15	Manufacturer's warranty	42

1 INTRODUCTION

Dear friends, thank you for choosing our products!

Non-contact thermometer for children **AGU NC8/NC8C/NC8A** has a wide measurement range, which allows for using it both for measuring body temperature and for measuring:

- milk surface temperature in a baby bottle;
- water surface temperature in a baby bathtub;
- ambient temperature.

Innovative technology allows for measuring temperature on the forehead in 1 second (less than 3 seconds for **NC8C/NC8A** models) at a distance up to 5 cm.

The unique design of the device with a built-in Heilmann infrared sensor provides accurate and reliable measurement results.

The temperature can be measured even in a sleeping child without disturbing.


The device automatically displays the last measurement result for 2 seconds after switching on.

5 paired beeps and red display backlight warn the user when the body temperature is equal to or above 37.5 °C.

Please read this User Manual carefully before using the thermometer.

2 LIST OF SYMBOLS


Symbol	Meaning
	The product meets the basic requirements of the Medical Devices Directive 93/42/EEC.
	WEEE (Waste Electrical and Electronic Equipment Directive). The symbol on the product or its package means that this product does not fall under the category of domestic waste. To avoid possible damage to the environment and human health, separate such wastes from others and dispose of them in accordance with accepted standards.
	The product complies with the basic requirements of TR CU 020/2011 «Electromagnetic compatibility of technical means».
	Warning/Attention.
	Read the User Manual before use.
	Type BF applied part.
	Manufacturer.
	Authorized manufacturer's representative in the European Union.

3 SCOPE OF APPLICATION

Non-contact thermometer for children **AGU NC8/NC8C/NC8A** is intended for domestic and clinical use. With the device you can measure body temperature on the forehead and any surfaces. This thermometer can be used both for children and adults. Do not use the device for purposes other than intended.

4 COMPLETE SET

- 1 Non-contact thermometer for children **AGU NC8/NC8C/NC8A** – 1 pc.
- 2 AAA battery – 2 pcs.
- 3 Storage case – 1 pc.
- 4 User Manual with Quick Guide – 1 pc.

5 BASIC SPECIFICATIONS

Type	Infrared digital non-contact electronic medical thermometer
Model	NC8/NC8C/NC8A
Power supply	3 V (2 AAA 1.5 V batteries)
Measurement range	Body temperature 32.0 °C ~ 42.5 °C (89.6 °F ~ 108.5 °F) Surface temperature 2–99 °C (35.6 °F ~ 199.0 °F)
Measuring accuracy	Object temperature: ±0.3 °C (±0.6 °F): 2.0–99.9 °C (35.6 °F ~ 199.0 °F) Body temperature: ±0.1 °C (±0.2 °F): 34.0 °C ~ 42.0 °C (93.2 °F ~ 107.6 °F) Out of range: ±0.3 °C (±0.5 °F)
Dimensions, mm	NC8 149x37x35; NC8C/NC8A 153x37x35
Weight (without batteries), g	NC8 ~ 67.2; NC8C ~ 65.6; NC8A ~ 67.0
Housing material	ABS plastic
Measurement time, sec	Approx. 1 second (around 3 for NC8C/NC8A models)
Display	Liquid crystal display, 4 characters with special icons, color backlight
Auto power-off	In 1 minute
Operating conditions	16–40 °C, humidity 15–65%
Memory	Automatic display of the last measured temperature. 9 memory cells for NC8A model

Indication when measuring body temperature

Temperature	Indication
$T < 32.0 \text{ °C (89.6 °F)}$	«L»
$32.0 \text{ °C (89.6 °F)} \leq T \leq 42.5 \text{ °C (108.5 °F)}$	Shows temperature value
$T > 42.5 \text{ °C (108.5 °F)}$	«H»


Acoustic signals

The thermometer is switched on and ready for use	1 short beep
Measurement completed	5 beeps in pairs if the value is above 37.5 °C
	2 short beeps if the value is below 35 °C
	1 long beep when measurement is complete
System error or malfunction	3 short beeps

6 DEVICE DESCRIPTION


Non-contact thermometer for children **AGU NC8** features:

- non-contact measurement method;
- temperature measurement within 1 second;
- high measurement accuracy;
- color and sound indication.


Non-contact thermometer for children **AGU NC8C** features:

- non-contact measurement method;
- temperature measurement within 2,5 seconds;
- high measurement accuracy;
- color and sound indication;
- distance indicator light.


Non-contact thermometer for children **AGU NC8A** features:

- non-contact measurement method;
- temperature measurement within 2,5 seconds;
- high measurement accuracy;
- color and sound indication;
- 9 memory cells;
- wireless communication module for synchronization with **AGU** app.

Battery indicator

7

GUIDELINES FOR SAFE OPERATION

Follow the instructions precisely to ensure reliable and long-term operation of the device.


- Any modification of the device is prohibited.
- Do not expose the thermometer to any mechanical shock or use if damaged.
- Do not immerse the device in water or other liquids.
- Do not use the device near strong electromagnetic sources, such as wireless communication sources or cellular phones.
- Do not allow children to play with the thermometer.
- Do not expose the thermometer to direct sunlight, extreme temperatures, moisture or contamination.
- The environmental temperature in the apartment may differ in various rooms. Before measuring temperature, the user and the thermometer should stay in the same room (premises) with constant environmental conditions for at least 15 minutes.
- Before starting measurement, please remove any dirt, hair or sweat from the measuring area.
- Do not measure temperature during or immediately after breastfeeding.
- The user should not take food or drink immediately before or during temperature measurement.
- Do not pocket the thermometer before use as it can be heated.
- Every time before starting measurement, please make sure that the infrared sensor lens is not contaminated. Clean the measuring sensor lens if it is contaminated.

**WARNING**

During early stages of illness with a rapid rise in temperature to high values, the effect of «white hyperthermia» can be observed –

a condition which can be characterized by the peripheral vessel constriction, when the skin becomes pale and remains cold. In such cases, temperature measurements on the forehead should not be carried out, because the skin temperature will be low.

8 BEFORE YOU START

- 1 Remove the battery compartment cover.
- 2 Install the batteries according to the instructions on the case.
- 3 Close the battery compartment cover.

9 OPERATION PRINCIPLE AND PROCEDURE

BATTERY INDICATOR

Charged battery



Battery low

Discharged battery


MEASURING BODY TEMPERATURE WITH AGU NC8/NC8C

- 1 To measure body temperature, turn on the device by pressing the  button once.
- 2 The display will show the result of the previous measurement (if the measurement was made earlier) and the thermometer will switch to stand-by mode. 
- 3 Bring the thermometer to the center of the forehead so that there is approximately 5 cm between the infrared sensor and the skin.

MEASURING WITH NC8


- 4 Press the **START** button once and release.

The display will show the body temperature measurement result. The thermometer will signal the end of the measurement with a


MEASURING WITH NC8C/NC8A

- 4 Press the **START** button to activate the distance indicator light. Focus the indicator light beam on one point as shown below.


The correct distance is reached.

beep and the display backlight will turn on.

The thermometer is too close.

The thermometer is too far away.

The measurement will be taken. The display will show the body temperature measurement result. The thermometer will signal the end of the measurement with a beep and the display backlight will turn on.


- 5** When temperature reaches 37.4 °C (or below), green backlight will turn on.

When temperature reaches 37.5 °C (or higher), red backlight will turn on.

OBJECT TEMPERATURE MEASUREMENT

- 1 To measure object temperature, turn on the device by pressing the  button once.
- 2 The display will show the result of the previous measurement (if the measurement was made earlier).
- 3 Bring the thermometer to the measuring area so that there is a distance of up to 5 cm between the sensor and the surface.
- 4 Press the **START** button 2 times. The thermometer will switch to object measurement mode and measure the temperature.
- 5 The display will show the object temperature measurement result. The thermometer will signal the end of the measurement with a beep and the display backlight will turn on.

SETTING THE DATE AND TIME FOR AGU NC8A

When using the thermometer for the first time and after each battery change, the date and time during the basic setting is displayed in the format 1-1 00: 00.

Perform step-by-step setting of the device:

- Press and hold the **ON/OFF** button for 5 sec.
- Select the year (20 to 35) by pressing the **START** button.
- Press the **ON/OFF** button again.
- Select the month (1 to 12) by pressing the **START** button.
- Press the **ON/OFF** button again.
- Select the day (from 1 to 31) by pressing the **START** button.
- Press the **ON/OFF** button again.
- Select the hour (0 to 23) by pressing the **START** button.

- Press the **ON/OFF** button again.
- Select the minutes (0 to 59) by pressing the **START** button.
- Press the **ON/OFF** button again to complete the setting.

PLAYING PREVIOUS DATA FOR AGU NC8A

Up to 9 data measurements are stored in the memory of the device.

- 1 Press the **START** button.
- 2 The display will show the memory cell number, then the data.
- 3 Press the **START** button to change the memory cells

When writing new data measurements (after 9), the old ones are cleared automatically.

TURNING ON/OFF THE LIGHT INDICATION AND SOUND NOTIFICATION

Go to the settings mode to enable/disable the light indication and sound notification.

- 1 To do this, hold the **START** button for 3 seconds with the thermometer turned off.
- 2 You will see the configurable parameter icon  and status **ON** (beep on) or **OFF** (beep off) on the display.
- 3 Press the **START** button to enable/disable the sound notification.
- 4 Wait for 3 seconds and the thermometer will switch to the light indication edit mode.
- 5 You will see the configurable parameter icon  and status **ON** (backlight and distance indicator light* are enabled) or **OFF** (backlight and distance indicator light* are disabled) on the display.
- 6 Press the **START** button to enable/disable the light indication.

- 7 After you have selected the desired mode, wait 4 seconds. The thermometer will save the selected parameters and turn off automatically. The next time you turn on the thermometer, it will work according to the preset parameters.

* For NC8C/NC8A models.

SWITCHING BETWEEN °C AND °F

To switch between measurement modes:

- 1 Turn off the thermometer.
- 2 Press and hold the **START** button for 8 seconds. The thermometer will switch to the measurement unit setting mode.
- 3 Press the **START** button to select the desired measurement mode in °C or °F.
- 4 After you have selected the desired mode, wait 4 seconds. The thermometer will save the selected parameters and turn off automatically. The next time you switch on the thermometer, it will work according to the preset parameters.

AGU NC8A OPERATION VIA MOBILE APPLICATION

AGU mobile application provides opportunity to:

- save the temperature measurement history;
- keep records of temperature data for different users;
- take notes on symptoms and general state;
- keep records of medications taken.

WARNING

The object temperature measurement history cannot be stored in the application.

Application installation: Android 5.0 or higher, iOS 10.0 or higher.


Option a: open the PlayMarket or AppStore, enter «AGU» in the search box and download the application.

Option b: for direct link to download the application, scan this QR code.

- 1 Download and run the application.
- 2 Register or log in through social media site.
- 3 Select device to connect.
- 4 Enter the required user data.
- 5 Switch on the device and press «Connect» button.

- 6 Follow the application instructions.

10 DEVICE CLEANING

WARNING

Do not use chemically active detergents to clean the device.

- 1 Clean the device with a soft, dry cloth. If necessary, the cloth can be slightly dampened.
- 2 If any dirt gets on the infrared sensor, the readings may differ from the actual ones, as this affects the measurement (usually the readings are too low). To eliminate this prob-

lem, it is enough to wipe the infrared sensor with a soft lint-free cloth (or cotton swab) moistened with disinfectant solution. After the cleaning procedure, the infrared sensor lens surface should be crystal clear. As the lens surface gets cool when the disinfectant solution evaporates, the thermometer should be left for several minutes at room temperature before use after cleaning.

- 3 Leave the device to air dry at room temperature.
- 4 Store the thermometer in a dry, clean place in accordance with the operating conditions.

11 TROUBLESHOOTING

Description	Possible cause	Solution
Temperature data are incorrect	Sensor is too close/too far from the measured surface	The distance from the measured surface to the sensor should be around 5 cm
Device won't turn on	Batteries are discharged	Replace the batteries
	Device is damaged	Contact the service center
	The ambient temperature is less than +16 °C	The ambient temperature should be in the range from +16 °C to +40 °C
	The ambient temperature exceeds +40 °C	
Incorrect (low) readings	IR sensor lens is dirty	Clean in accordance with section 10 of this Manual

12 STORAGE, TRANSPORTATION AND OPERATION RULES

- The device should be transported at temperatures ranging from $-10\text{ }^{\circ}\text{C}$ to $+40\text{ }^{\circ}\text{C}$ and a relative humidity of up to 75%.
- The device should be operated at temperatures ranging from $+16\text{ }^{\circ}\text{C}$ to $+40\text{ }^{\circ}\text{C}$ and a relative humidity of up to 65%.
- The device should be stored at temperatures ranging from $-20\text{ }^{\circ}\text{C}$ to $+50\text{ }^{\circ}\text{C}$ and a relative humidity of up to 80%.
- Do not expose the device to thermal shock.

WARNING

After transportation or storage at low temperatures, it is necessary to keep the device at room temperature for at least 2 hours before switching on. Remove the batteries if the thermometer is not going to be used for a long time.

13 DISPOSAL

The symbol on the product or its package means that this product does not fall under the category of domestic waste.

- Proper disposal of the device will prevent adverse environmental and human health effects.
- In order to protect the environment, the device must not be disposed of together with domestic (household) waste. Disposal shall be provided in accordance with local regulations.
- The device must be disposed of in accordance with the EU Directive 2012/19/EU WEEE (Waste Electrical and Electronic Equipment).
- If you have any questions, please contact the local public utility responsible for waste disposal.

14 CERTIFICATION

Non-contact thermometer for children **AGU NC8/NC8C/NC8A** conforms to the provisions of the EU Medical Device Directive 93/42/EEC.

15 MANUFACTURER'S WARRANTY

Warranty information is specified in the warranty card. The warranty period starts from the date of sale, providing that all operation conditions specified in this User Manual are strictly observed. The warranty is valid only on presentation of the sales receipt and the warranty card completed by the authorized dealer. The warranty does not apply to consumables, constituent parts and components, such as battery, package, etc.

Warranty and free maintenance service is not provided in case of:

- Use of the device in violation of the operating instructions.
- Damage caused by intentional or erroneous actions of the consumer due to improper or negligent handling.
- Evidence of mechanical impact, dents, cracks, chips, etc. on the thermometer housing, evidence of opening the device, disassembly, unauthorized repair, ingress of moisture, exposure to aggressive agents or any other unauthorized alterations of the device, and in other cases of violation of storage, cleaning, transportation and operation rules specified in the User Manual.
- Ingress of oils, dust, insects, liquids and other foreign objects inside the device.

Follow the instructions precisely to ensure reliable and long-term operation of the device. In case of abnormal operation of the device, please contact the service center or retailer.

The manufacturer reserves the right to make alterations in the design of the device.

WARNING

Do not attempt to repair the device yourself – this will void the warranty.

The thermometer alterations without the manufacturer's permission will void the warranty.

For repair and maintenance, please contact a specialized service center (for more information, please visit www.agu-baby.com).

Date of manufacture and importer information are specified on the individual package.

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AGU[®] is the registered trademark by Montex Swiss AG,
MedNet EC-REP GmbH,
Shenzhen Combei Technology Co., Ltd.