Ray

# SMART THERMOMETER

NON-CONTACT  
AGU SNC10

de

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## CONTENTS

1	Introduction.....	27
2	List of symbols.....	28
3	Scope of application .....	29
4	Complete set.....	29
5	Basic specifications .....	30
6	Device description.....	31
7	Guidelines for safe operation .....	33
8	Before you start .....	34
9	Operation principle and procedure .....	35
10	Cleaning .....	40
11	Troubleshooting .....	41
12	Storage, transportation and operation rules .....	42
13	Disposal.....	42
14	Certification .....	43
15	Manufacturer's warranty .....	46

## 1 INTRODUCTION

Dear friends, thank you for choosing our products!

Smart non-contact thermometer **AGU SNC10** has a wide measurement range, which allows you to use it both for measuring body temperature and for measuring:

- milk surface temperature in a baby bottle;
- water surface temperature in a baby bathtub;
- ambient temperature.

Innovative infrared sensor technology allows for measuring temperature on the forehead in 1 second.

The unique design of the device with a built-in infrared sensor provides accurate and reliable measurement results.

It is easy and simple to measure temperature with smart non-contact thermometer **AGU SNC10**.

Having a distance sensor, **AGU SNC10** starts measuring automatically as soon as the correct distance between the thermometer and the measured surface is reached.

The temperature can be measured in a sleeping child without disturbing him.

This thermometer has the ability to synchronize with a mobile application. This allows you to store measurement results and add notes to measurements.

Please read this user manual carefully before using the device.

## 2 LIST OF SYMBOLS

Symbol	Meaning
	CE marking with Notified Body identification number. Indicates compliance with the European Medical Device Directive 93/42/EEC.
	The product complies with the basic requirements of TR CU 020/2011 «Electromagnetic compatibility of technical means».
	WEEE (Waste Electrical and Electronic Equipment Directive). The symbol on the product or its package means that this product does not fall under the category of domestic waste. To avoid possible damage to the environment and human health, separate such wastes from others and dispose of them in accordance with accepted standards.
	Special disposal of batteries. The symbol on the product or its package indicates that the batteries are not classified as household waste. To avoid possible damage to the environment and human health, separate such wastes from others and dispose of them in accordance with accepted standards.
	Type BF device (body floating) – Type B device with an insulated F-Type applied part. Type of protection of applied part against electric shock.
<b>IP22</b>	Ingress protection rating. Leading digit (protection against ingress of solid foreign objects): 2 – protection against

ingress of solid objects more than 12 mm in size; fingers or other objects with a maximum length of 80 mm, or solid objects. Second digit (protection against ingress of foreign liquids): 2 – protection against vertically dripping water at an angle of up to 15° from the vertical (equipment in normal position).

Warning/Attention.

Read the user manual before use.

Manufacturer.

Authorized manufacturer's representative in the European Union.


### 3 SCOPE OF APPLICATION

Smart non-contact thermometer **AGU SNC10** is intended for domestic and clinical use. With the device you can measure body temperature on the forehead and any surfaces. This thermometer can be used both for children and adults. Do not use the device for purposes other than intended.

### 4 COMPLETE SET

- 1 Smart non-contact thermometer **AGU SNC10** – 1 pc.
- 2 1.5 V AAA battery – 2 pcs.
- 3 User manual with quick guide – 1 pc.

2 X 1.5 V AAA

## 5 BASIC SPECIFICATIONS

Type	Infrared non-contact thermometer
Model	AGU SNC10
Power supply	3 V (2 x 1.5 V AAA batteries)
Measurement range	Body temperature: 34.0 °C ~ 43 °C (93.2 °F ~ 109.4 °F) Object temperature: 2 °C – 93.2 °C (35.6 °F ~ 199.9 °F)
Measuring accuracy	Object temperature: $T \leq 40 \text{ °C} \pm 2 \text{ °C} (\pm 0.4 \text{ °F})$ ; $T > 40 \text{ °C} \pm 5\%$ Body temperature: $\pm 0.2 \text{ °C} (\pm 0.4 \text{ °F})$ : 34.0 °C ~ 42.0 °C (93.2 °F ~ 107.6 °F) Out of range: $\pm 0.3 \text{ °C} (\pm 0.5 \text{ °F})$
Dimensions, mm	~147 (L) x ~38 (W) x ~48 (H)
Weight (with batteries), g	Approx. 100
Housing material	ABS plastic, acryl
Measurement time	Around 1 sec
Display	Segment display with special icons
Auto power-off	30 seconds
Operating conditions	Temperature: 15 °C ~35 °C (59 °F ~ 95 °F) Humidity: up to 95% (non-condensing)
Memory	9 measurements

### Indication when measuring body temperature

Temperature	Indication
$T < 34.0 \text{ °C} (93.2 \text{ °F})$	«L»
$34.0 \text{ °C} (93.2 \text{ °F}) \leq T \leq 43.0 \text{ °C} (109.4 \text{ °F})$	Shows temperature value
$T > 43.0 \text{ °C} (109.4 \text{ °F})$	«H»

## Acoustic signals


The thermometer is switched on and ready for use	1 short beep
Measurement completed	5 beeps if the value is above 37.5 °C
	1 long beep if the value is below 37.5 °C
System error or malfunction	3 short beeps

## 6 DEVICE DESCRIPTION

Smart non-contact thermometer **AGU SNC10** features:


- non-contact measurement method;
- temperature measurement within 1 second;
- high measurement accuracy;
- Illumination of the measurement area;
- sound indication.

Wireless synchronization with AGU app.

- 1 IR sensor.
- 2 Battery compartment cover.
- 3 Segment display.
- 4 **START** button.
- 5 **ON\OFF** button.
- 6 Illumination of the measurement area.
- 7 Distance sensor.

## DESCRIPTION OF THE DISPLAY ELEMENTS

- 1 Measurement result in Fahrenheit.
- 2 Measurement result in Celsius.
- 3 Low battery level.
- 4 Forehead measurement mode.
- 5 Measurement in process.
- 6 Object temperature measurement mode.
- 7 Sound notification enabled.
- 8 Measurement result.

## **7** GUIDELINES FOR SAFE OPERATION

Follow the instructions precisely to ensure reliable and long-term operation of the device.

- Any modification of the device is prohibited.
- Do not expose the thermometer to mechanical shock or use if damaged.
- Do not immerse the device in water or other liquids.
- Do not use the device near strong electromagnetic sources, such as wireless communication sources or cellular phones.
- Do not allow children to play with the thermometer.
- Do not expose the thermometer to direct sunlight, extreme temperatures, moisture or contamination.
- The environmental temperature in the apartment may differ in various rooms. Before measuring temperature, the user and the thermometer should stay in the same room (premises) with constant environmental conditions for at least 15 minutes.
- Before starting measurement, please remove any dirt, hair or sweat from the measuring area.
- Do not measure temperature during or immediately after breastfeeding.
- The user should not take food or drink immediately before or during temperature measurement.
- Do not pocket the thermometer before use as it can be heated.
- Every time before starting measurement, please make sure that the infrared sensor lens is not contaminated. Clean the measuring sensor lens if it is contaminated.


## WARNING

During early stages of illness with a rapid rise in temperature to high values, the effect of «white hyperthermia» can be observed – a condition which can be characterized by the peripheral vessel constriction, when the skin becomes pale and remains cold. In such cases, temperature measurements on the forehead should not be carried out, because the skin temperature will be low.

## 8 BEFORE YOU START

- 1 Remove the device from the package.
- 2 Remove the battery compartment cover.
- 3 Insert new 2 x 1.5 V AAA batteries into the battery compartment, observing polarity. Always use only new batteries.
- 4 Close the battery compartment cover.

Dispose of used batteries in accordance with applicable legal regulations. It is forbidden to dispose of batteries in household waste.

## 9 DEVICE SETTINGS

### DATE AND TIME SETTINGS

Before the first use of the thermometer and after each replacement of batteries, date and time are displayed in the format 1-1 00:00.

Perform a step-by-step setup of the device:

- Press and hold the **ON\OFF** button for 5 seconds.
- Set the year (from 20 to 35) by pressing the **START** button.  
Press the **ON\OFF** button again.
- Set the month (from 1 to 12) by pressing the **START** button.  
Press the **ON\OFF** button again.
- Set the day (from 1 to 31) by pressing the **START** button.  
Press the **ON\OFF** button again.
- Set the hour (from 0 to 23) by pressing the **START** button.  
Press the **ON\OFF** button again.
- Set the minutes (from 0 to 59) by pressing the **START** button.  
Press the **ON\OFF** button again to complete the setup.

### REPRODUCTION OF THE PREVIOUS READINGS


You can recall up to 9 measurement values.

- 1 Press the **START** button.
- 2 The display will show the memory location number, then the value.
- 3 Press the **START** button to switch memory locations.

When recording new measurement values (after the 9th value), old values are deleted automatically.

## SOUND ON/OFF

- 1 Press and hold the **START** button for 5 seconds.
- 2 The display will show the beeper icon and its status.


Press the **START** button to turn the sound notification **ON** or **OFF**

- 3 Do not press any buttons on the device for 3 seconds. The thermometer will turn off after saving the settings.

## SWITCHING BETWEEN THE UNITS OF MEASUREMENT

- 1 Press and hold the **START** button for 9 seconds.
- 2 The display will show °C or °F.

- 3 Press the **START** button to switch between the units of measurement.
- 4 Do not press any buttons on the device for 3 seconds. The thermometer will turn off after saving the settings.

## OPERATION PRINCIPLE AND PROCEDURE

There are 2 options for the device operation:

**Option 1:** without mobile application.

**Option 2:** via mobile application.

## DEVICE OPERATION WITHOUT MOBILE APPLICATION

Smart non-contact thermometer **AGU SNC10** measures the body temperature on the forehead as well as the temperature of objects.

### MEASURING BODY TEMPERATURE ON THE FOREHEAD


- 1 Remove any dirt, hair or sweat before placing the temperature sensor in the measurement area.
- 2 Turn on the thermometer by pressing the **ON\OFF** button.

- 3 The thermometer will go into standby mode.

The thermometer is ready for measurement.

- 4 Bring the thermometer to the center of the forehead at a distance of no more than 3 cm and press the **START** button. The arrow on the display will flash during measurement. During the measurement, the illumination of the measurement area is also activated.


The temperature will be measured if the distance to the surface is no more than 3 cm.

- 5 When the measurement is complete, you will hear a beep and the display will show the measured body temperature.

## MEASURING OBJECT TEMPERATURE

- 1 Turn on the thermometer by pressing the **ON\OFF** button.
- 2 The thermometer will go into standby mode.  
The thermometer is ready for measurement.
- 3 Bring the thermometer to the measurement surface at a distance of no more than 3 cm and press the **START** button twice.

- 4 You will hear a beep and the display will show the measured object temperature.

## DEVICE OPERATION VIA MOBILE APPLICATION

**AGU** mobile application provides opportunity to:

- save the temperature measurement history;
- save the temperature measurement history for different users;
- take notes on symptoms and general state;
- take notes on medicines taken.

### **WARNING**

The object temperature measurement history is not stored in the application.

**Application installation:** Android 6.0 or higher, iOS 13.0 or higher.


**Option a:** open the PlayMarket or AppStore, enter «AGU» in the search box and download the application.

**Option b:** for direct link to download the application, scan this QR code.

- 1 Download and run the application.
- 2 Register or log in through social media site.
- 3 Select device to connect.
- 4 Enter the required user data.
- 5 Switch on the device and press «Connect» button.

- 6 Follow the application instructions.


## 10 DEVICE CLEANING

### WARNING

Do not use chemically active detergents to clean the device.



- 1 Clean the device with a soft, dry cloth. If necessary, the cloth can be slightly dampened.
- 2 If any dirt gets on the infrared sensor, the readings may differ from the actual ones, as this affects the measurement (usually the readings are too low). To eliminate this problem, it is enough to wipe the infrared sensor with a soft lint-free cloth (or cotton swab) moistened with disinfectant solution. After the cleaning procedure, the infrared sensor lens surface should be crystal clear. As the lens surface gets cool when the disinfectant solution evaporates, the thermometer should be left for several minutes at room temperature before use after cleaning.

- 3 Leave the device to air dry at room temperature.
- 4 Store the thermometer in a dry, clean place in accordance with the operating conditions.

## 11 TROUBLESHOOTING

If a problem occurs during operation, first of all, refer to the list of possible failures given below.

Indication	Cause	Solution
	The measured temperature is higher: in body temperature measurement mode – $T > 43.0\text{ °C}$ (109.4 °F)	Use the thermometer only for the specified temperature ranges: If necessary, clean the sensor tip. In the event of a repeated error message, contact your retailer or Customer Service
	The measured temperature is higher: in object temperature measurement mode – $T > 93.2\text{ °C}$ (or 199.9 °F)	
	The measured temperature is lower: in body temperature measurement mode – $T < 34.0\text{ °C}$ (93.2 °F)	
	The measured temperature is lower: in object temperature measurement mode – $T < 2.0\text{ °C}$ (35.6 °F)	
	1. Operating temperature is out of range between 15 °C – 35 °C (59 °F – 95 °F). 2. Too frequent temperature measurements	

## 12 STORAGE, TRANSPORTATION AND OPERATION RULES

- The device should be stored at temperatures between  $-25\text{ }^{\circ}\text{C}$  and  $+55\text{ }^{\circ}\text{C}$  ( $-13\text{ }^{\circ}\text{F}$  ...  $+131\text{ }^{\circ}\text{F}$ ) with a maximum relative humidity of 95% (non-condensing).
- The device should be operated at temperatures between  $+15\text{ }^{\circ}\text{C}$  and  $+35\text{ }^{\circ}\text{C}$  ( $59\text{ }^{\circ}\text{F}$  –  $95\text{ }^{\circ}\text{F}$ ) with a relative humidity of up to 95% (non-condensing).
- The device should be transported at temperatures between  $-25\text{ }^{\circ}\text{C}$  and  $+55\text{ }^{\circ}\text{C}$  ( $-13\text{ }^{\circ}\text{F}$  ...  $+131\text{ }^{\circ}\text{F}$ ) with a maximum relative humidity of 85% (non-condensing).
- Do not expose the device to thermal shock.

### WARNING

After transportation or storage at low temperatures, it is necessary to keep the device at room temperature for at least 2 hours before switching on.

## 13 DISPOSAL

The symbol on the product or its package means that this product does not fall under the category of domestic waste.

- With proper disposal of the product, you will help to prevent possible negative effects of the device on the environment and human health.
- For the protection of the environment, the device should not be disposed of with domestic (household) waste. Disposal must be carried out in accordance with local regulations.

- Device must be disposed of in accordance with the EU Directive 2012/19/EU WEEE (Waste Electrical and Electronic Equipment).

If you have any questions, please contact the local public utility responsible for waste disposal.

Dispose of used, fully discharged batteries in special containers, take to special waste collection centers or electronics stores. Users are legally obliged to ensure the disposal of batteries.

## 14 CERTIFICATION

This product complies with the provisions of Medical Devices Directive (MDD) 93/42/EEC.

EMC Tables		
<p><b>AGU SNC10</b> is intended for use in the electromagnetic environment specified below. The customer or the user of <b>AGU SNC10</b> must make sure that it is used in such an environment</p>		
<p>Guidance and manufacturer's declaration – Electromagnetic emissions</p>		
Phenomenon	Professional healthcare facility environment a)	HOME HEALTHCARE ENVIRONMENT a)
Conducted and radiated RF EMISSIONS	a)	CISPR 11 Group 1 Class B
Harmonic distortion	Not applicable	

Voltage fluctuations and flickering	Not applicable		
<p>a) The equipment is suitable for use in Home Health Environments and Professional Health Care Environments limited to patient rooms and respiratory treatment facilities in hospital or clinics. The more restrictive acceptance limits of Group 1 Class B (CISPR 11) have been considered and applied. The equipment is suitable for use in the mentioned environments when directly connected to the Public Mains Network.</p> <p>b) The test is not applicable in this environment unless the ME EQUIPMENT and ME SYSTEM used will be connected to the PUBLIC MAINS NETWORK and the power input is otherwise within the scope of the Basic EMC standard</p>			
Guidance and manufacturer's declaration – Electromagnetic immunity – Enclosure port			
Phenomenon	Basic EMC standard or test method	Immunity test levels	
		Professional health-care facility environment	HOME HEALTHCARE ENVIRONMENT
ELECTROSTATIC DISCHARGE	IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	

Radiated RF EM fields	IEC 61000-4-3	a)	10 V/m b) 80MHz –2.7 GHz 80% AM at 1kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	<p>COMPLIANT</p> <p><b>NOTE:</b> Further information about distances to be maintained between portable and mobile RF communications equipment (transmitters) and the <b>AGU SNC10</b> can be requested from supplier using the contact information provided in this manual. However, it is advisable to keep the electro-mechanical aerosol equipment at an adequate distance of, at least, 0.5 m from mobile phones or other RF communications transmitters to minimize possible interference</p>	

RATED power frequency magnetic fields	IEC 61000-4-8	30 A/m c 50 Hz or 60 Hz
<p>a) The equipment is suitable for use in Home Health Environments and Professional Health Care Environments limited to patient rooms and respiratory treatment facilities in hospital or clinics. The more restrictive IMMUNITY acceptance limits have been considered and applied.</p> <p>b) Before modulation is applied.</p> <p>c) This test level assumes a minimum distance of at least 15 cm between the ME EQUIPMENT or ME SYSTEM and sources of power frequency magnetic fields</p>		

## 15 MANUFACTURER'S WARRANTY

Warranty information is specified in the warranty card. The warranty period starts from the date of sale, providing that all operation conditions specified in this User Manual are strictly observed. The warranty is valid only on presentation of the sales receipt and the warranty card completed by the authorized dealer. The warranty does not apply to consumables, constituent parts and components, such as battery, packaging, etc.

Free warranty service is not provided in case of:

- use of the device in violation of the operating instructions;
- damage caused by intentional or erroneous actions of the consumer due to improper or negligent handling; evidence of mechanical impact, dents, cracks, chips, etc. on the thermometer housing, evidence of opening the device, disassembly, unauthorized repair, ingress of moisture, exposure to aggressive agents or any other unauthorized alterations of the device, and in other cases of viola-

tion of storage, cleaning, transportation and operation rules specified in the User Manual;

- ingress of oils, dust, insects, liquids and other foreign objects inside the device.

Follow the instructions precisely to ensure reliable and long-term operation of the device.

In case of abnormal operation of the device, please contact your retailer or Customer Service.

The manufacturer reserves the right to make alterations in the design of the device.

**WARNING**

Do not attempt to repair the device yourself – this will void the warranty.

The thermometer alterations without the manufacturer's permission will void the warranty.

For repair and maintenance, please contact a specialized after-sales service (see [agu-baby.com](http://agu-baby.com)).

Date of manufacture and importer information are specified on the individual package.

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Medical Device Safety Service GmbH (MDSS),
AViTA Corporation, 9F, No. 78, Sec. 1, Kwang-Fu Rd.,