Tomchi

# NEBULIZER

MESH

AGU N8

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EN USER MANUAL

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**CONTENT**

1	Introduction .....	97
2	List of symbols .....	98
3	Scope of application .....	100
4	Completeness .....	101
5	Specifications .....	103
6	Components .....	107
7	Safe operation instructions .....	108
8	Preparation for work .....	110
9	Operating principles and procedures .....	112
10	Troubleshooting.....	126
11	Storage, transportation and operation conditions .....	128
12	Disposal .....	129
13	Certification .....	130
14	Manufacturer's warranty .....	139

# 1

## INTRODUCTION

Dear friends, thank you for choosing our products!

Mesh nebulizer **AGU N8** is a medical device. Please follow the physician's prescriptions to choose the right type, dosage and regime. Spray characteristics vary depending on the drug used. The nebulization rate may vary depending on the medicine.

### WARNING

As with any electronic device, this product may fail due to a power outage, battery discharge or mechanical failure.

Please read this instruction manual carefully before using the device.

## 2

## LIST OF SYMBOLS

## Symbols Meaning

- The device complies with CE requirements for nebulizers, inhalers, separators and converters.
- The product meets the requirements of Directive 93/42/EEC concerning medical equipment.

WEEE (EC Directive on Waste Electrical and Electronic Equipment). This marking means that the device should not be disposed of together with household waste. In order to prevent a possible damage to the environment and harm to human health, such waste should be segregated from other types of waste and disposed of in accordance with the relevant accepted standards.

The unit meets the requirements of Customs Union Technical Regulation 020/2011 on Electromagnetic Compatibility of Technical Products.

Classification:

- Equipment with internal power source
- Product of type BF
- IP22
- Not suitable for use in the presence of a flammable anesthetic mixture with oxygen or nitrous oxide
- Continuous operation mode.

IP22

The degree of protection against insertion. The first figure (protection against insertion of foreign solid objects): 2 - protection against insertion of solid objects, the size of more than 12 mm; fingers or other objects with a length of not more than 80 mm, or solid objects. The second figure (protection against foreign liquids): 2 - protection against falling drops, falling objects from above at an angle to the vertical no more than 15 ° (the device in the normal position).

Authorized Representative of the manufacturer in the European Union.

In order to avoid the incorrect operation of the device caused by electromagnetic interference between electrical and electronic equipment, do not use the device near mobile phones or microwave ovens.

Warning/Note/Attention.

Read the instruction before use.

Protect from direct sunlight.

Manufacturer.

### 3

## SCOPE OF APPLICATION

### Intended use

Certified medical staff such as doctors, nurses and therapists, medical personnel, or patients under the supervision of qualified medical staff.

The user should also understand the operating principle of **AGU N8** and read the instruction manual.

### Intended user

Adults and children with asthma, chronic obstructive pulmonary diseases (COPD), such as emphysema and chronic bronchitis, or other respiratory diseases that are characterized by airway obstruction.

### Usage environment

This product is intended for use in medical institutions, such as hospitals, clinics and medical rooms, in living rooms or outdoors under a roof.

### Operation period

The life of the device is 24 months, if the device is used to nebulize the saline solution 3 times a day for 10 minutes each time at room temperature (23 °C). The length of the operating period may vary depending on the environment of use.

The nebulizer is a medical device. Please follow the physician's prescriptions to choose the right type, dosage and regime.

### Nebulization rate of the device

Nebulization rate  $\geq 0.35$  ml/min (may vary depending on drugs).


- Thank you for purchasing this product. To ensure the safe and correct use of the device, be sure to read the user's manual.
- Please keep this instruction manual in an accessible place for future reference.
- This device is a unit of individual use.
- Don't allow several patients to use the same device without prior replacement or disinfection of a mouthpiece or a mask.
- The tip and mask can be reused after disinfection.

## 4


**COMPLETENESS**

The kit includes the following components. If you find that there are no components, immediately contact the authorized seller from whom you purchased the device.

Main unit

Protective cover

Medication chamber

Batteries


AA 1.5Vx2 (optional)


Tip/mouthpiece


Carrying bag


Micro - USB  
cable (optional)

Inhalation mask (S)

Inhalation mask (M)


\* Liquid condenser

\* The liquid condenser is designed to collect a drug in a volume of 4 ml.

## Main components

## 5

## SPECIFICATIONS

<b>Product name</b>	<b>Mesh Nebulizer</b>
Model	AGU N8
Nebulization method	Active Vibrating Mesh (Active Vibrating Membrane)
Dimensions, LxWxH, mm	~ 65×54×113
Weight, g	~ 96 (without batteries)
Power source	-3V DC («AA» 1.5V alkaline battery x2) -Micro - USB socket 5V DC 1A
Vibration frequency, kHz	~ 113
Nebulization rate, ml/min	≥ 0,35*
Particle size, µm	MMAD 3,6
Recommended chamber filling volume, ml	~ 8 maximum ~ 0.5 minimum
Battery life	Battery life up to 3.0 hours with continuous use. 10 days if used daily for 20 minutes (10 minutes 2 times a day). (Use 2 «AA» (LR6) alkaline batteries)
Device service life	The service life specified below is applicable in case the device is used for spraying the physiological solution for 10 minutes 3 times per day at room temperature (23 °C). The specified service life of the device may depend on the environment in which it is used. Main unit – 24 months. Medication chamber – 12 months
Operating position	Vertically. Up to 40° inclination within 15 seconds
Lack of medicine indicator	The lack of drugs indicator flashes slowly (blue) (about 1 Hz) for 15 seconds and then the power turns off automatically
Warranty	The warranty* for the main unit makes 24 months, while for the medication chamber – 6 months as from the date of acquisition. The warranty* is valid subject to availability of a warranty card filled out by a duly authorized representative confirming the date of sale and also of the relevant receipt. The warranty does not cover the consumable parts, components and constituent elements, namely: protective cap, battery, tip / nozzle, carrying bag, micro-USB cable, face masks (S, M), condensate pot, package. *For more details about the warranty service conditions refer to Section 14 "Manufacturer's Warranty".
Operation	10~40 °C (50~104 °F), 15~93% relative humidity
Storage	-20~70 °C (-4~158 °F), ≤ 93% relative humidity
Accessories	Protective cover, mouthpiece, alkaline batteries (optional), carry bag, user manual, inhalation mask (S), inhalation mask (M)

\* When tested with saline, at a normal temperature of 23°C and a constant current of DC 5V, 0.5 A (when using a Micro - USB adapter).

The nebulizer generates a high-frequency sound and automatically turns off if the medicine does not enter the membrane of the medication chamber for more than 15 seconds (the time varies for different types of solutions) or if the drug is running out. It helps to prevent the membrane damage.

## SPECIFICATIONS: (TECHNICAL DATA)


### NOTE

The test results reflect cascade impact measurement of the sprayed particles.

Productivity may vary depending on different types of medicines, such as suspensions or solutions with high viscosity. For more information, see the list of medicines' data.

MMAD = Mass median aerodynamic diameter.

The particle size corresponds to EN 13544-1.

Picture 4

Cumulative curve based on the data in Table 3

Table 2 summarizes the test results, and Table 3 shows the mass fraction of the medicine at each stage. The graph of cumulative size distribution is shown in Picture 4.

**Table 2**  
**The summarized test results for the indicated device (Mesh nebulizer, model AGU N8)**

Test	Specifications	
Determination of particle size	Total mass of emitted particles (M1, $\mu\text{g}$ )	$7392.59 \pm 122.47$
	The mass of particles collected by the cascade impactor (M2, $\mu\text{g}$ )	$297.90 \pm 15.18$
	The mass of particles collected by an external filter (M3, $\mu\text{g}$ )	$1746.5 \pm 31.84$
	Inactive volume in the medication chamber ( $\mu\text{g}$ )	$42607.41 \pm 122.47$
	Recovery (%) <sup>1</sup>	$27.66 \pm 0.71$
	MMAD ( $\mu\text{m}$ )	$1.53 \pm 0.04$
	Cum. val. MMAD for four tests (%)	2.34
	Geometric Standard Deviation (GSD)	$1.7 \pm 0.02$
	Cum. val. GSD for four tests (%)	0.93
	Inhaled fraction (% , 0.5-5 $\mu\text{m}$ )	$94.69 \pm 0.29$
	Inhaled mass ( $\mu\text{g}$ , 0.5-5 $\mu\text{m}$ )	$282.10 \pm 15.2$
	Procedure time (min)	1
	Fraction of large particles (%) (> 4.7 $\mu\text{m}$ )	$5.81 \pm 0.32$
	Fraction of small particles (%) (> 4.7 $\mu\text{m}$ )	$94.19 \pm 0.32$
Fraction of ultra-small particles (%) (< 1 $\mu\text{m}$ )	$17.25 \pm 1.44$	
Output power and speed of aerosol delivery	Aerosol delivered capacity (ml) <sup>2</sup>	$1.74 \pm 0.05$
	Aerosol delivery rate (ml / min)	$0.27 \pm 0.01$

<sup>1</sup> Recovery is calculated by the formula  $(M2 + M3) / M1$ .

<sup>2</sup> Continue the procedure until the medication chamber becomes empty or until the formation of steam stops.

**Table 3**  
**Mass fraction of the medicine at each stage of the impactor influence**


Stage No.	Effective cut diameter	Cumulative fraction (%)					Average value per minute at each stage ( $\mu\text{m}/\text{in}$ )	Average amount of fraction at each stage (%)
		Test 1	Test 2	Test 3	Test 4	Average value		
1	21.30	100.00	100.00	100.00	100.00	100.00	6.75	2.27
2	14.80	97.77	97.49	97.63	98.03	97.73	4.50	1.51
3	9.80	96.50	95.69	96.27	96.38	96.21	9.50	3.19
4	6.00	93.63	92.10	93.22	93.09	93.01	24.75	8.31
5	3.50	85.04	84.55	84.75	84.54	84.72	85.05	28.55
6	1.60	58.31	56.82	54.92	54.61	56.16	136.55	45.84
7	0.93	11.46	11.57	9.49	8.88	10.35	18.00	6.04
8	0.53	4.46	4.74	3.73	4.28	4.30	8.50	2.85
After filtering	<53	1.27	1.51	1.02	1.97	1.44	4.30	1.44
Total							297.90	

## 6


## COMPONENTS

The kit includes the following components. If you find that there are no components, immediately contact the seller from whom you purchased the device.

Main unit

Protective cover



Medication chamber

Alkaline batteries  
AA 1.5Vx2 (optional)


Tip/mouthpiece


Carrying bag

Micro - USB  
cable (optional)

Inhalation mask (S)

Inhalation mask (M)

\* Liquid condenser

\* The liquid condenser is designed to collect a drug in a volume of 4 ml.

# 7

## SAFE OPERATION INSTRUCTIONS

In order to ensure safe and correct use of the device, read the manual carefully before use.

### WARNING

- Please follow the physician's prescriptions to choose the right type, dosage, and regime.
- Do not use any solutions in the medication chamber that are not prescribed by your physician. This is a device for individual use. Do not allow several patients to use the same device without first replacing or disinfecting the mouthpiece or mask.
- If you use a nebulizer after purchasing it or you have not used it for a long time, please clean the sprayer parts. (See page 116-119).
- Rinse the medication chamber, silicone condenser and tip with distilled water after each use. Immediately dry the parts and store them in a clean place. (See page 116-119).
- Before the first use, the inhalation mask should be rinsed with distilled water and dried.
- The tip and mask can be reused after disinfection (see how to disinfect on page 118-119).

### USAGE PRECAUTIONS

- If the device does not turn off automatically after out of drug and produces a high-frequency sound, press the «I/O» button to turn off the power immediately to avoid the membrane damage. For malfunctions detection, go to page (124).
- Please clean the nebulizer parts carefully after each use. Otherwise, it may not work.
- Water is not allowed to use. If you fill the medication chamber with water, the nebulizer cannot be turned on. Distilled water can be used to clean the medication chamber as part of the «cleaning process».
- Do not use cotton buds or any foreign objects to contact the membrane of the medication chamber. It may damage the membrane.

- Do not drop the device. Avoid strong impact on the nebulizer.
- Do not mix different types of batteries.
- Do not store or transport the nebulizer if there is medicine or water left in it.
- Do not allow the contact water with the main unit of the nebulizer.
- Keep the device away from children. Children should only use the device under adult supervision.

### **Product characteristics**

- Compact and light.
- Low energy consumption and a little medicine residue.


The nebulizer can function normally for a short time after turning the device at any angle. If the nebulizer is turned so that the medicine does not enter the membrane, the nebulizer can spray properly for about 15 seconds. (The spraying time varies depending on the specific types of medicine).

## 8 PREPARATION FOR WORK

### HOW TO ASSEMBLE THE NEBULIZER

#### 1. Attach the medication chamber to the main unit:

- a) Install the liquid condenser (1.1).
- b) Insert the medication chamber until hearing a click sound into the main unit (1.2).
- c) Instal the cover into the main unit (1.3).


Make sure that the medication chamber is installed correctly; otherwise, it may lead to a bad connection, and the nebulizer may not work properly.

Please keep the electrodes of the main unit and the medication chamber clean; otherwise, the nebulizer may not work properly.

## 2. Attach the mouthpiece or inhalation mask.

- Tightly attach the tip to the main unit.
- Please clean and dry the parts of the nebulizer before use. (See page -116).


You can either use batteries or a DC 5V power supply as a power source for this product.


### How to install batteries

#### 1. Open the battery compartment and install 2 AA alkaline batteries.

a) Open the battery compartment cover.

b) Insert the batteries so that the polarity is oriented correctly as shown in the picture.

c) Close the battery compartment cover.

## Battery life and replacement

- Up to 3 hours of continuous use. Use 2 AA alkaline batteries (LR6) for 6 days, when using 3 times a day for 10 minutes.
- If the power indicator (orange color) blinks 2 times per second, then the battery is almost discharged. However, the nebulizer can still be used for about 30 minutes.
- If the power indicator lights constantly (orange color), the nebulizer does not work because of low battery. Replace alkaline batteries immediately.


### **WARNING**

Do not use different types of batteries.

The battery life may vary depending on the batteries used.

## 2. Press the on/off button (I/O) (1).


- Power indicator:
  - a) The device operates normally if the indicator permanently lights green.
  - b) The battery is partially discharged if the indicator blinks in orange 2 Hz.
  - c) The battery is fully discharged if the indicator lights constantly orange.

## HOW TO FILL THE CHAMBER (CONTAINER) WITH DRUGS

### 1. Remove the medication chamber from the main unit:

A small flange  
to take cover away

- take off the cover by pulling a small flange to take a cover away. (1.1);
- press the PUSH button to remove the medication chamber from the main unit (1.2).

### **!** WARNING

Please, in order to avoid damaging the nebulizer, be sure to press the PUSH button before you start removing the medication chamber.

Please, in order to prevent a rupture of a membrane, do not touch it by finger or other objects.


## 2. Fill the medication chamber with drugs.

- Fill in the medicine as shown in the picture. Recommended filling capacity: minimum 0.5 ml/not more than 8 ml.
- Close the lid of the medication chamber.

### **WARNING**


Please, in order to avoid leakage of the drug from the chamber, make sure that the lid of the chamber is securely closed.

Always fill the chamber with the drug when it is removed from the main unit.

## 3. Reattach the medication chamber to the main unit.

Attach the chamber to the main unit.

### **WARNING**

Make sure that the medication chamber is installed correctly, as incorrect installation may cause malfunction.

Make sure that the medication chamber is installed correctly; otherwise, problems in the device operation may arise.

Please keep the electrodes of the main unit and the medication chamber clean; otherwise, problems in the device operation may arise.

Clean the medication chamber after each use.


## HOW TO USE A NEBULIZER

Before assembling the device, you can pour 0.9% sodium chloride solution (table salt) into the medication chamber before use, and then press the «I/O» button. If the nebulizer does not spray the liquid, refer to part «Troubleshooting» on page (123).

### 1. Turn on the power.

a) Press the «I/O» button, the power indicator should constantly light (green). The backlight of the medication chamber should constantly light (blue) (1.1).

b) Medication indicator with a blue light blue light flash with 1Hz if no medication (empty & depleted) (1.2).

### **!** WARNING


If you turn on the device with an empty medication chamber, the power indicator will turn green for a short time before it goes out. After that, the illumination of the medication chamber will start blinking blue for about 15 seconds, then the device will automatically turn off.

During operation, if the medication is finished in the medication chamber, the power indicator will turn off, and the illumination of the medication chamber will start blinking blue at a frequency of 1 time per sec. for about 15 seconds, then the illumination will turn off and the device will automatically turn off.

**Maximum continuous operating time of the nebulizer:** 30 minutes.

## 2. Inhalation.

Holding the nebulizer securely in your hand to start the inhalation.

### WARNING

If the device does not detect any medicine in the container, it will automatically turn off.

If the device does not turn off automatically when there are no medicine, press the «I/O» button to turn off the power immediately, to avoid damaging the membrane. To determinate malfunctions, please go to page 124.

During inhalation, the nebulizer can be kept at any angle, however, make sure that the medicine is in contact with the membrane, otherwise, the device will turn off automatically after 15 seconds.


At the end stage of the medication, it is recommended to tilt the device slightly (front side) towards yourself. This will allow the remaining medicine to reach the membrane for spraying.

Do not shake the nebulizer strongly during use; otherwise, the device may automatically turn off.

Provide careful supervision when the child uses the nebulizer.

### 3. Turn off the power.

- The nebulizer automatically turns off after 15 seconds after out of drugs.
- If you want to stop the inhalation, press the «I/O» button.
- Press the «I/O» button.

### DEVICE CLEANING AFTER USE

Be sure to clean the medicine cup with distilled water after each use before transportation or storage.

#### **WARNING**

#### 1. Remove any medicine residues.

- Open the lid of the medicine cup and remove any medicine residues.
- Use citric acid solution to clean the medicine cup\*. Method of preparation of citric acid solution: pour 100 ml of drinking water into a clean container or glass. Add half a teaspoon of citric acid (approx. 2.5 - 3.0 g) to 100 ml of drinking water and mix thoroughly. Pour the prepared citric acid solution into the medicine cup to a level of 4 ml and close the medicine cup lid.

\*Use crystalline or powdered citric acid to prepare the solution.

- Press and hold the ON/OFF button for 5 seconds. The device will turn on and switch to cleaning mode, and the LED indicator of the medicine cup will start flashing blue at a rate of 3 times per second, the network indicator will start flashing orange against a background of constant green light. Spray a citric acid solution for two minutes to remove the medicine residues. If necessary, repeat the above procedure until the medicine residues are completely removed.


 If the citric acid solution in the medicine cup runs out and the device produces a high-frequency sound, press the ON/OFF button to turn off the device. Otherwise, the membrane may wear out.

Please remove the medicine residues after each use. Otherwise, the medicine cup membrane may become clogged.

- Rinse the medicine cup with distilled water after cleaning with citric acid solution. To do this, pour the distilled water into the medicine cup to a level of 6 ml, close the lid and shake the cup at least 5 times. Then open the lid and pour out the water. Rinse the outer part of the membrane with distilled water to remove citric acid residues. Repeat the cleaning procedure with distilled water at least 3 times.

## 2. Disassemble the nebulizer.

- Remove the inhalation mask or mouthpiece from the nebulizer (2.1).
- Remove the cover (2.2).
- Remove the medication chamber from the nebulizer (2.3).


## WARNING

Before the first use, the inhalation mask should be rinsed with water and dried. Some of the components of the device can be re-used repeatedly after disinfection (see how to disinfect on page 119).


## 3. Rinse the parts with enough distilled water.

Rinse the medication chamber, mouthpiece, inhalation mask and liquid condenser with water.

#### 4. Dry the cleaned parts.

After cleaning, the parts should be thoroughly dried with clean gauze cloth.

#### **!** WARNING

Please do not use cotton cloth or fabric from other materials to clean the membrane; otherwise, dust or fabric fiber may remain on the membrane, leading to malfunction of the nebulizer.


Please do not allow any cotton buds or foreign objects to contact with the membrane of the medication chamber.

#### 5. Wipe the main unit with a clean gauze cloth.

- Gently wipe the main unit by a wet gauze. Then, dry the device by using a new, clean, dry gauze.


#### **!** WARNING

Please clear the electrodes on the main unit and the medication chamber. It will ensure normal electrical conductivity and, therefore, normal spraying.

6. Attach the medication chamber to the main unit and close the cover of the main unit. Store all the details in a clean place.


## 7. Disinfection.


### **⚠ WARNING**

In the case of infectious diseases, the device should be disinfected every time after use and all parts should be sterilized. You can use alcohol or other common disinfectant sprays based on alcohol for disinfection. Make sure that there are no remaining disinfectant residues on the parts to ensure safe inhalation on next use.

Disinfect the tip/mouthpiece by boiling at 100 °C for 3 minutes.

Please be careful, do not boil the medication chamber and mask, as they may be damaged.

Carefully dry the parts

## REPLACEMENT OF THE MEDICATION CHAMBER

The service life of a medication chamber is 12 months under normal conditions (when used three times a day for a total use of 30 minutes). However, spraying efficiency may begin to deteriorate in less than 12 months depending on how you use it or when using certain types of medications. If the nebulizer does not work correctly or the spraying speed is significantly reduced after cleaning, you should replace the medication chamber with a new one. (If you want to purchase a medication chamber, contact the authorized seller from whom you purchased the product or the nearest distributor).


### 1. Remove the medication chamber from the nebulizer.

Open/disassemble the main unit. Press the PUSH button on the back of the main unit and pull the container towards the front of the main unit.

#### WARNING

In order to avoid the device damage, make sure that you press the PUSH button before removing the container from the main unit.

In order to avoid the membrane rupture, do not press on it or poke it with your finger or other objects.


a) Take off the cover by pulling a small flange to take a cover away. (1.1).

b) Press the PUSH button to remove the medication chamber from the main unit (1.2).

## 2. Installation of the medication chamber.

Install the medication chamber correctly as shown in the picture.

- Install the liquid condenser (2.1).
- Attach the chamber to the main unit (2.2).
- Close the cover of the main unit (2.3).

Please make sure that the medication chamber is installed correctly, otherwise, due to a bad connection, there may be problems with the device.

Please keep the contacts of the main unit and medication chamber clean, otherwise, the problems in the operation of the device may arise.


Please rinse the medication chamber before use.

## HOW TO TRANSPORT A NEBULIZER

In order to transport the nebulizer correctly, it is necessary to disassemble the parts and perform the actions described below. Then place the accessories in the carrying bag.


### 1. Disassemble the nebulizer.

Please remove the mouthpiece or inhalation mask as shown in the picture.

### 2. Install the protective cover.

Please install the protective cover as shown in the picture. Ensure the panel is well installed. It will protect the nebulizer from possible damage during transport.


### 3. Place the main unit and parts in the carrying bag.

**⚠ WARNING**

Please do not transport a nebulizer that still contains medicines or water. Liquids can leak, as well as stain or damage the device.

Do not store the nebulizer in a place with high temperature, humidity or direct sunlight.

## 10 TROUBLESHOOTING

In case of malfunctions in the device operation, refer to the following table for its elimination.

<b>Malfunctions</b>	<b>Possible reasons</b>	<b>Troubleshooting</b>
Extremely low nebulizing speed	The medication chamber is not well installed	Install the medication chamber correctly and start nebulizing. (See page 108)
	The drugs do not enter the membrane for more than 15 seconds	Place the nebulizer at such an right angle that the drugs get on flow to the membrane
	The membrane of the medication chamber is clogged	Clean the medication chamber. If this does not solve the problem, replace a new chamber
	The electrodes on the medication chamber are littered with medicine or water	Clean the electrodes and start nebulizing again. (See page 118)
	The electrodes on the nebulizer and medication chamber are contaminated	Clean the electrodes and restart the power. (See page 118)
After turning on, the power indicator lights up for 15 seconds and then turns off immediately	The medication chamber is not well installed	Install the medication chamber correctly and start nebulizing. (See page 101)
	The medication chamber is empty	Fill in the medication chamber. (See page 113)
	The medicine does not enter the membrane	Place the nebulizer at such angle that the can get on the membrane
	The electrodes on the nebulizer and medication chamber are contaminated	Clean the electrodes and restart the power. (See page 118)
The power indicator is off and the nebulizer does not work	The batteries are installed incorrectly	Reinstall the batteries, observing the polarity and restart the device. (See page 110)
	Low battery	Replace the batteries and restart the device. (See page 110)
	The micro - USB cable is not properly connected to the nebulizer	Connect it correctly and restart the device. (See page 110)

The power indicator is on, and the nebulizer does not work	If the low battery indicator lights, it means that the batteries are low or the power is insufficient to operate the device	Replace the batteries and reinstall the device. (See page 110)
	The membrane of the medication chamber is damaged	Replace the medication chamber with a new one and fill in the drugs. (See page 113)
	The contacts electrodes on the medication chamber are contaminated with medicine or water	Clean the contacts electrodes from the remnants of drugs or water and restart the device. (See page 118)
	The electrodes on the nebulizer or medication chamber are contaminated	Clean the electrodes and restart the device. (See page 118)
	The membrane of the medication chamber is strongly clogged	If it still not work normally after cleaning, replace the medication chamber with a new one.
The nebulizer turns off during the use	The medication chamber is not installed up to the stop	Install the medication chamber correctly and restart the device. (See page 113)
	The micro - USB cable is not properly connected to the nebulizer	Connect it correctly and restart the device. (See page 109)
	Out of drugs	Add the drugs to the medication chamber. (See page 113)
	The drugs do not enter the membrane for more than 15 seconds	Place the nebulizer at right angle that the drugs flow to the membrane
The nebulizing process is unstable	The nebulizer was shaken when using	Hold the nebulizer securely in your hand. (See page 115)
The nebulizer does not turn off automatically when out of the drugs	The medication chamber is broken	Replace the medication chamber with a new one and fill in the drugs. (See page 113)
	Some solutions for inhalation can lead to the formation of a large amount of foam in the medication chamber	Clean the foam and restart the device. (See page 113)
	The contacts electrodes on the medication chamber are contaminated with medicine or water	Clean the contacts electrodes and restart the device. (See page 118)
	The electrodes on the nebulizer or medication chamber are contaminated	Clean the electrodes and restart the device. (See page 118)
Leakage of drugs from the medication chamber	Damage to the medication chamber or run out of the silicone gasket	Replace the medication chamber with a new one and fill in the drugs. (See page 113)

If your nebulizer still does not work properly (go to [agu-baby.com](http://agu-baby.com)).

## **11 STORAGE, TRANSPORTATION AND OPERATION CONDITIONS**

The device should be transported at a temperature of  $-20^{\circ}\text{C}$  to  $+70^{\circ}\text{C}$  with a maximum relative humidity of 75% or less.

The device can be operated at a temperature of  $+10^{\circ}\text{C}$  to  $+40^{\circ}\text{C}$  with a relative humidity of 15 ~ 93%.

Store the device at a temperature of  $-20^{\circ}\text{C}$  to  $+70^{\circ}\text{C}$  with a maximum relative humidity of not more than 93%.

Do not subject the device to sudden temperature fluctuations.

After transport or storage at low temperatures, the device should be kept at room temperature for at least 2 hours before switching on.

## 12

## DISPOSAL

The symbol on the device or on its packaging indicates that the product is not classified as household waste.

- With the proper disposal of the product, you will help to prevent possible negative effects of the device on the environment and human health.
- The device should not be disposed of with domestic (household) waste for the environment protection. The product should be disposed of in accordance with local regulations.
- The device should be disposed of in accordance with the Directive EU 2012/19/EU WEEE (Waste Electrical and Electronic Equipment) for waste electrical and electronic equipment.

If you have any questions, please contact the local utility responsible for waste disposal.

## 13 CERTIFICATION

### Electromagnetic compatibility

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.

### Electromagnetic emanation: instructions and statements of the manufacturer.

The device is intended for use in the electromagnetic environment described below.

Buyers and users of the device must ensure that it is used in an environment that meets the guidelines.

Emanation test	Compliance	Statements on the electromagnetic environment
HF emanation CISPR 11	Group 1	High frequencies are used only to provide internal functions of the device. That is why the level of HF emanation is very low, and it is un-likely that it will interfere with near-by electronic equipment
HF emanation CISPR 11	Class B	The device is suitable for use at home, as well as inside any premises that are directly connected to a standard electrical network de-signed to supply electricity to the housing stock
Emanation of harmonic components of IEC 61000-3-2	Class A	
Fluctuations/pulsations/ voltages IEC 61000-3-3	Compliance	

### Electromagnetic emanation: instructions and statements of the manufacturer.

The device is intended for use in the electromagnetic environment described below.

Buyers and users of the device must ensure that it is used in an environment that meets the guidelines.

<b>Immunity test</b>	<b>IEC 60601 Test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment-guidance</b>
Electrostatic discharge IEC 61000-4-2	6 kV contact 8 kV air	6 kV contact 8 kV air	The floor must be wooden, concrete or covered with ceramic tiles. If the floor is covered with synthetic materials, the relative humidity should be at least 30%
Fast electrical transients/ flashes IEC 61000-4-4	2 kV for electrical power lines 1 kV for input/output lines	2 kV for electrical power lines 1 kV for input/output lines	The quality of electricity in the power network must meet the standards for power supply systems for industrial or medical environments
Flashes IEC 61000-4-5	1 kV between the lines 2 kV between the line (s) and the ground	1 kV between the lines 2 kV between the line (s) and the ground	The quality of electricity in the power network must meet the standards for power supply systems for industrial or medical environments
Interruptions and voltage variations on the input power lines IEC 61000-4-11	<5% UT (UT drop >95%) during the half cycle  40% UT (UT drop 60%) during 5 cycles  70% UT (UT drop 30%) during 25 cycles  <5% UT (UT drop >95%) during 5 seconds	<5% UT (UT drop >95%) during the half cycle  40% UT (UT drop 60%) during 5 cycles  70% UT (UT drop 30%) during 25 cycles  <5% UT (UT drop >95%) during 5 seconds	The quality of electricity in the power network must meet the standards for power supply systems for industrial or medical environments. If the user needs the continuous operation of the device in the presence of voltage interruptions, it is recommended to use batteries or an uninterruptible power supply to power the device
Magnet fields of the power frequency (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	The quality of electricity in the power network must meet the standards for power supply systems for industrial or medical environments
NOTE: UT is the voltage in the supply network before the test begins			

# ELECTROMAGNETIC COMPATIBILITY: INSTRUCTIONS AND STATEMENTS OF THE MANUFACTURER - FOR A DEVICE THAT IS NOT A LIFE SUPPORT DEVICE

## Electromagnetic compatibility: instructions and statements of the manufacturer

The device is intended for use in the electromagnetic environment described below.

Buyers and users of the device must ensure that it is used in an environment that meets the guidelines.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment-guidance
HF conductivity IEC 61000-4-6 HF emanation IEC 61000-4-3	3 V (avg.) 150 kHz – 80 MHz 3 V/m 80MHz – 2,5 GHz	3 V (avg.)  3 V/m	Portable and mobile radio communication equipment should not be used closer to any part of the device (including cables) than the recommended minimum distance calculated with the appropriate equation for the frequency of the transmitter of such equipment  $d = 1,2\sqrt{P}$ $d = 1,2\sqrt{P}$ 80 MHz – 800 MHz $d = 2,3\sqrt{P}$ 800MHz – 2,5GHz  Where P is the maximum transmitter power in watts (W) specified by the manufacturer and d is the recommended minimum distance in meters (m)  The measured values of the electromagnetic force produced by the stationary radio transmitter <sup>a</sup> should be less than the corresponding allowable values for all frequency bands <sup>b</sup> . There may be interference near equipment with the following symbol 

### NOTE 1

At 80 MHz and 800 MHz, the values for the high-frequency range are used.

### NOTE 2

These instructions may be applicable not in all situations. The propagation of electromagnetic waves is affected by their absorption and reflection by structures, objects, and people.

a. The value of the electromagnetic force generated by stationary transmitters, such as base stations of mobile phones and mobile radio systems, amateur radio transmitters, AM and FM radio transmitters and television transmitters, cannot be exactly determined theoretically. To assess the possible impact of such a transmitter, it is necessary to carry out electromagnetic measurements on site. If the measured values of the electromagnetic field at the site of use of the **AGU N8** exceed the corresponding permissible level, it is necessary to check the correct operation of the **AGU N8**. If the device does not operate normally, additional measures may be taken, such as changing the orientation of the **AGU N8** or its location.

b. In the frequency range 150 kHz - 80 MHz, the field force values should be less than 3 V/m.

**The recommended minimum distance between portable and mobile radio equipment and the device - for a device, that is not a LIFE SUPPORT DEVICE!**

**The recommended minimum distance between portable and mobile radio equipment and the device.**

The device is intended for use in an electromagnetic environment in which the level of emitted high-frequency interference is monitored. Users or purchasers of the device can help prevent the occurrence of electromagnetic interference by observing the minimum permissible distances between the portable or mobile high-frequency communication equipment (transmitter) and the device, as recommended below, depending on the maximum emanation power of the communication equipment.

Nominal maximum emanation power of the transmitter (W)	The minimum distance, depending on the frequency of the transmitter (M)		
	150 kHz – 80 MHz $d = 1,2\sqrt{P}$	80 MHz - 800 MHz $d = 1,2\sqrt{P}$	800 MHz – 2,5 GHz $d = 2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters with nominal maximum emanation power not specified above, the recommended minimum distance  $d$  in meters (m) can be estimated using the appropriate transmitter frequency of the equation, where  $P$  is the maximum transmitter power in watts (W) according to the transmitter manufacturer's data.

## NOTE 1

At 80 MHz and 800 MHz, the values for the high-frequency range are used.

## NOTE 2

These instructions may be applicable not in all situations. The propagation of electromagnetic waves is affected by their absorption and reflection by structures, objects, and people.

### Appendix A: EMC information

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

#### Guidance and Manufacturer's declaration – electromagnetic emissions

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

## Guidance and Manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	6kV contact 8 kV air	6kV contact 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %
Electrical fast transient/burst IEC 61000-4-4	2kVfor power supply lines 1kVfor input/output lines	2kVfor power supply lines 1kVfor input/output lines	Mains power quality should be that of a typical commercial or hospital environment
Surge IEC 61000-4-5	1 kV line(s) to line(s) 2kV line(s) to earth	1 kV line(s) to line(s) 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment
interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle  40 % UT (60 % dip in UT) for 5 cycles  70 % UT (30 % dip in UT) for 25 cycles  <5 % UT (>95 % dip in UT) for 5 sec	<5 % UT (>95 % dip in UT) for 0,5 cycle  40 % UT (60 % dip in UT) for 5 cycles  70 % UT (30 % dip in UT) for 25 cycles  <5 % UT (>95 % dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment
NOTE UT is the AC mains voltage prior to application of the test level			

## GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY – FOR DEVICE THAT IS NOT LIFE SUPPORTING

### Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter  <b>Recommended separation distance</b>  $d = 1,2$ $d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800 MHz to 2,5 GHz  Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	

#### NOTE 1

At 80 MHz and 800 MHz, the higher frequency range applies.

#### NOTE 2

These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

### Recommended separation distances between portable and mobile RF communications equipment and the device – for device that is not LIFE SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the device			
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1,2$	80 MHz to 800 MHz $d = 1,2$	800 MHz to 2,5 GHz $d = 2,3$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

## **NOTE 1**

At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

## **NOTE 2**

These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## 14 MANUFACTURER'S WARRANTY

The warranty for the main unit makes 24 months, for the medication chamber – 6 months as from the date of acquisition.

The warranty is valid subject to availability of a warranty card filled out by a duly authorized representative confirming the date of sale and also of the relevant receipt. The warranty does not cover the consumable parts, components and constituent elements, namely: protective cap, battery, tip / nozzle, carrying bag, micro-USB cable, face masks (S, M), condensate pot, package.

### **The warranty and free maintenance service are not provided in case:**

- the device is used with violation of the requirements set by the operation manual;
- damage is caused by consumer's willful or wrong actions as a result of improper or negligent use;
- traces of mechanical effect (dents, cracks, chips, etc.), tampering, dismantling, traces of repairs outside the authorized technical service center, traces of ingress of moisture into the case or effect of corrosives, or of any other outside interference with the device design are observed on the device case, and also in other cases of violation by the consumer of the rules of storage, cleaning, transportation and maintenance of the device stipulated by the operation manual;
- penetration of oil, dust, insects, liquids (not intended for the use with devices) and other foreign matters inside the device.

### **The warranty does not cover the defects (faults) of the device caused by the following reasons:**

- normal wear and tear of component parts with a limited service life;
- damage caused to accessories and nozzles which do not form an integral part of the device (protective cap, battery, tip / nozzle, carrying bag, micro - USB cable, face masks (S, M), condensate pot, package);
- deposition of sediment on the medication chamber (grid / membrane) irrespective of the used inhalation solution;
- use of low-quality accessories and nozzles, or the accessories and nozzles which reached the end of their service life;
- force majeure circumstances (accident, fire, flood, electrical power failures, etc.).

## **WARNING**

Strictly follow the instructions in order to ensure reliable and long-term operation of the device.

Apply to the sales assistance in case the device does not work properly.

For the repair and maintenance apply to a specialized service department.

The manufacturer reserves a right to change the design of the device.

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MedNet GmbH, Borkstrasse 10, 48163 Münster, Germany.