Giraffe

THERMOMETER

NON-CONTACT
FOR CHILDREN

AGU NC8

DE

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FR

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EN

USER MANUAL
Version 1.0

RU

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CONTENT

| | | |
|----|---|----|
| 1 | Introduction..... | 33 |
| 2 | List of symbols..... | 34 |
| 3 | Scope of application | 35 |
| 4 | Completeness | 35 |
| 5 | Specifications..... | 36 |
| 6 | Device description..... | 37 |
| 7 | Safe operation instructions..... | 38 |
| 8 | Preparation for work..... | 39 |
| 9 | Operating principles and procedures..... | 39 |
| 10 | Device cleaning..... | 42 |
| 11 | Troubleshooting..... | 43 |
| 12 | Storage, transportation and operation terms..... | 44 |
| 13 | Disposal..... | 44 |
| 14 | Certification | 45 |
| 15 | Manufacturer's warranty..... | 45 |

1 INTRODUCTION

Dear friends, thank you for choosing our products!

Non-contact thermometer for children **AGU NC8** has a wide measurement range, which allows it to be used as a thermometer for measuring body temperature, as well as:

- milk surface temperature in a baby bottle;
- water surface temperature in a baby bath;
- ambient temperature.

The latest technology with the use of an infrared sensor allows measuring temperature on forehead in 1 second.

The unique design of the device with built-in infrared sensor provides accurate and reliable measurement results.

Non-contact thermometer for children **AGU NC8** allows conducting temperature measurements easily and simply.

The temperature can be measured even in a sleeping child without disturbing.

The measurement is fast, which is very convenient when taking a child's temperature.

The device automatically displays the last measurement for 2 seconds after switching on.

5 paired beeps and a red backlight of LCD warns the user that the body temperature is equal to or higher than 37.5°C.

Please read the instruction manual carefully before using the device.

2 LIST OF SYMBOLS


| Symbol | Meaning |
|---|---|
|  | The product meets the requirements of Directive 93/42/EEC concerning medical equipment. |
|  | WEEE (Waste Electrical and Electronic Equipment Directive). The symbol on the product or on its packaging indicates that this product is not classified as household waste. Please, separate such waste from other waste types and dispose it according to accepted standards, in order to avoid possible harm to environment and human health. |
|  | The unit meets the requirements of Customs Union Technical Regulation 020/2011 on Electromagnetic Compatibility of Technical Products. |
|  | Warning/Caution. |
|  | Read the instruction manual before use. |
|  | BF type product. |
|  | Manufacturer. |
|  | Authorized representative of the manufacturer in the European Union. |

3 SCOPE OF APPLICATION

Non-contact thermometer for children **AGU NC8** allows measuring body temperature on forehead and any surfaces. The device is designed for home use. Do not use the device for other purposes. The device is designed to measure temperature, both in children and adults.

4 COMPLETENESS

- 1 Non-contact thermometer for children **AGU NC8** – 1 pc.
- 2 AAA battery – 2 pcs.
- 3 Storage pouch – 1 pc.
- 4 Instruction manual – 1 pc.

5 SPECIFICATIONS

| | |
|-------------------------------|--|
| Type | Infrared digital non-contact electronic medical thermometer |
| Model | NC8 |
| Power source | 3 V (2 batteries AAA x 1.5 V) |
| Measurement range | Body temperature 32.0°C ~ 42.5°C (89.6°F ~ 108.5°F) Surface temperature 2–99°C (35.6°F ~ 199.0°F) |
| Measurement accuracy | Object temperature: ±0.3°C (±0.6°F): 2.0–99.9°C (35.6°F ~ 199.0°F) Body temperature: ±0.1°C (±0.2°F): 34.0°C ~ 42.0°C (93.2°F ~ 107.6°F) Out of range: ±0.3°C (±0.5°F) |
| Dimensions, mm | 149x36x35 |
| Weight (without batteries), g | Approx. 68.5 |
| Housing material | ABS plastic |
| Measurement time | 1 second |
| Display | Liquid crystal display, 4 digits with special icons, color illumination |
| Auto power off | In 1 minute |
| Operation terms | 16–40°C, humidity 15–65% |
| Memory | Automatic display the last measured temperature |

Indication of body temperature measurement

| Temperature | Indication |
|--|------------------|
| $T < 32.0^{\circ}\text{C}$ (89.6°F) | «L» |
| 32.0°C (89.6°F) $\leq T \leq 42.5^{\circ}\text{C}$ (108.5°F) | Show temperature |
| $T > 42.5^{\circ}\text{C}$ (108.5°F) | «H» |

Sound signals


| | |
|---|--------------|
| The device is switched on and ready for measurement | 1 short beep |
|---|--------------|

| | |
|-----------------------------|---|
| Measurement completion | 5 paired signals if the value is higher than 37,5°C |
| | 2 short beeps if the value is below 35°C |
| | 1 long beep when the measurement is completed |
| System error or malfunction | 3 short beeps |

6 DEVICE DESCRIPTION

The features of Non-contact thermometer for children **AGU NC8**:

- non-contact measurement method;
- temperature measurement in 1 second;
- high accuracy of measurements;
- color and sound indication.

7

SAFE OPERATION INSTRUCTIONS


The device must be used for its intended purpose, in accordance with this manual.

- Any modification of the device is prohibited.
- Do not expose the thermometer to mechanical influences and do not use it if damaged.
- Do not immerse the device in water or other liquids.
- Do not use the device near high electromagnetic radiation sources, for example, near wireless sources or cellular phones.
- Do not allow children to play with the thermometer.
- Do not expose the thermometer to direct sunlight, heat, moisture, or dirt.
- Replace the batteries when needed.
- The environmental temperature in the apartment may differ in various rooms. Before measuring temperature, the user and the thermometer should stay in the same room (premises) with constant environmental conditions for at least 15 minutes.
- Before starting measurement, please remove any dirt, hair or sweat from the measuring area.
- Do not measure temperature during or straight after nursing.
- The user should not take food or drink immediately prior or during temperature measurement.
- Do not pocket the thermometer before use as it can be heated.
- Every time before starting measurement, please make sure that the measuring sensor lens is not contaminated. In case of contamination, clean the measuring sensor lens.

During early stages of illness with a rapid rise in temperature to high values, the effect of «white hy-

perthermia» can be observed – a condition which can be characterized by the vessel constriction, when the skin becomes pale and remains cold. In such cases, temperature measurements on the forehead should not be carried out, because the temperature of the skin will be low.

8 PREPARATION FOR WORK


- 1 Open the battery cover.
- 2 Install the batteries according to the instructions on the case.
- 3 Close the battery cover.


9 OPERATING PRINCIPLES AND PROCEDURES

BATTERY INDICATION

BODY TEMPERATURE MEASUREMENT

- 1 To measure body temperature, switch on the device by pressing the button  once.
- 2 The result of the previous measurement will be shown on the display (when measuring previously). 
- 3 Put the thermometer to the center of the forehead, maintaining a distance of 3–4 cm between the sensor and the skin.
- 4 Press the **START** button once and release the button.

- 5 The result of the body temperature measurement appears on the device display. The thermometer will signal the measurement completion by a sound signal and by the display backlight.

At the temperature 37.4°C and below, the green backlight turns on.

At the temperature 37.5°C and above, the red backlight turns on.

OBJECT TEMPERATURE MEASUREMENT

- 1 To measure object temperature, switch on the device by pressing the button  once.
- 2 The result of the previous measurement will be shown on the display (when measuring previously).
- 3 Bring the thermometer to the measured surface so that there is a distance of 3–5 cm between the sensor and the surface.
- 4 Press the **START** button 2 times and release. The thermometer will be switched to the measurement mode and take temperature of the object.
- 5 The result of the object temperature measurement appears on the device display. The thermometer will signal the measurement completion by a sound signal and by the display backlight.

ACTIVATE OR INACTIVATE LIGHT AND SOUND NOTIFICATION

To activate the light and sound notification, it needs to enter the settings mode.

As the thermometer off, hold the **START** button for 3 seconds.

You will see the icon of the adjustable parameter on the device display  and the status **ON** (sound on) or **OFF** (sound off).

Press the **START** button to turn the sound notification on or off.

Wait for 3 seconds and the thermometer will switch to the light notification mode.

You will see the icon of the adjustable parameter on the device display  and the status **ON** (light on) or **OFF** (light off).

Press the **START** button to turn the light notification on or off.

Wait for 4 seconds after selecting the desired mode. The thermometer will save the parameters and turn off automatically. The next time you turn on, the thermometer will work according to the set parameters.

SWITCHING BETWEEN THE MEASUREMENT °C AND °F

You should turn off the thermometer to switch between the measurement modes.

Then press and hold the **START** button for 8 seconds. The thermometer will enter the measurement units setting mode.

Press the **START** button to select the desired measurement mode °C or °F.

Wait for 4 seconds after selecting the desired mode. The thermometer will save the parameters and turn off automatically. The next time you turn on the thermometer will work according to the set parameters.

10 DEVICE CLEANING

ATTENTION

Do not use chemically active detergents to clean the device.

- 1 Clean the device with a piece of dry cloth. If necessary, the cloth can be slightly moistened.
- 2 In case of contamination of the infrared thermometer sensor lens, the readings begin to change immediately as something interferes with the measurement (usually the readings are too low). To eliminate this problem, it is enough to wipe the sensor with a soft lint-free cloth (or cotton swab) moistened with disinfectant solution. Then it is neces-

sary to wipe the lens with a dry soft lint-free cloth. After the cleaning procedure, the lens surface should be crystal clear. As the lens surface gets cool when the disinfectant solution evaporates, the thermometer should be left for several minutes at room temperature before use after cleaning.

- 3 Leave the device to dry in air at room temperature.
- 4 Store the device in a dry, clean and dust-free place in accordance with the operating terms.

11 TROUBLESHOOTING

| Description | Possible reasons | Solution |
|---|---|--|
| Temperature data is not displayed properly | The sensor is too close/too far from the measured surface | The distance from the measured surface to the sensor should be 3–4 cm (for body temperature measuring) and 3–5 cm (for object temperature measuring) |
| The device can't switch on | Low battery | Replace the battery |
| | Device damage | Contact the service center |
|  | Ambient temperature is below +16°C | Ambient temperature should be in the range from +16°C to +40°C |
|  | Ambient temperature is higher than +40°C | |
| Wrong (low) readings | The IR sensor lens is dirty | Clean the lens according to item 10 of this instructions manual |

12 STORAGE, TRANSPORTATION AND OPERATION TERMS

- The device should be transported at -10°C to $+40^{\circ}\text{C}$ with a maximum relative humidity of 75% or less.
- The device can be operated at a temperature of $+16^{\circ}\text{C}$ to $+40^{\circ}\text{C}$ with a maximum relative humidity of not more than 65%.
- The device should be stored at a temperature of -20°C to $+50^{\circ}\text{C}$ with a maximum relative humidity of not more than 80%.
- Do not subject the device to sudden temperature fluctuations.

ATTENTION

After transport or storage at low temperatures, it is necessary to keep the device at room temperature for at least 2 hours before switching on.

It is recommended to remove the batteries if the device is not supposed to be used for a long time.

13 DISPOSAL

The symbol on the product or on its packaging indicates that this product is not classified as household waste.

- With proper disposal of the product, you will help to prevent possible negative effects of the device on the environment and human health.
- For the protection of the environment, the device should not be disposed of with domestic (household) waste. Disposal must be carried out in accordance with local regulations.

- The device must be disposed of in accordance with the EU Directive 2012/19/EU WEEE (Waste Electrical and Electronic Equipment) for waste electrical and electronic equipment.

If you have any questions, please contact the local utility responsible for waste disposal.

14 CERTIFICATION

Non-contact thermometer for children **AGU NC8** complies with Directive 93/42/EEC on medical equipment.

15 MANUFACTURER'S WARRANTY

The warranty period of this product is 12 months from the date of sale with strict observance of the operating terms in accordance with this instructions manual. The warranty does not apply to battery. Carefully follow the instructions to ensure reliable and long-term operation of the device.

If the device does not work as expected, contact the service center.

The manufacturer reserves the right to make changes in the design of the device.

Do not try to repair the device yourself; it will lead to loss of warranty.

Changes made to this thermometer without manufacturer's permission will void the warranty.

For repair and maintenance, please contact a specialized after-sales service (see agu-baby.com).

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AGU[®] is the registered trademark by Montex Swiss AG,
Shenzhen Combei Technology Co., Ltd.