Eaglet

THERMOMETER

INFRARED

AGU IHE3

DE

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EN

USER MANUAL

RU

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CONTENT

| | | |
|----|---|----|
| 1 | Introduction..... | 25 |
| 2 | List of symbols..... | 26 |
| 3 | Scope of application..... | 27 |
| 4 | Completeness..... | 27 |
| 5 | Basic specifications..... | 28 |
| 6 | Product description..... | 29 |
| 7 | Safe operation instructions..... | 30 |
| 8 | Preparation for work..... | 31 |
| 9 | Operation principle and procedure..... | 32 |
| 10 | Cleaning..... | 38 |
| 11 | Troubleshooting..... | 39 |
| 12 | Storage, transportation and operation requirements..... | 39 |
| 13 | Disposal..... | 40 |
| 14 | Certification..... | 40 |
| 15 | Manufacturer's warranty..... | 45 |

1 INTRODUCTION

Dear friends, thank you for choosing our products!

Infrared thermometer **AGU IHE3** is an advanced high-quality product that is compliant with international standards.


Utilizing infrared technology, Infrared thermometer **AGU IHE3** takes temperature in seconds by measuring heat generated by ear drum and surrounding tissues or skin surface of the forehead and surrounding tissues. Innovative infrared sensor technology ensures one-second temperature measurement. This device is suitable for people of all ages.

The thermometer can be also used for measuring ambient temperature.

Please read this manual carefully before using the device.

2 LIST OF SYMBOLS


| Symbol | Meaning |
|---|--|
|  | The CE marking with the identification number of the Notified Body. This denotes the compliance of European Medical Device Directive 93/42/EEC. |
|  | The product complies with the basic requirements of TR CU 020/2011 «Electromagnetic compatibility of technical means». |
|  | WEEE (Waste Electrical and Electronic Equipment Directive). The symbol on the product or its package means that this product does not fall under the category of domestic waste. To avoid possible damage to the environment and human health, separate such wastes from others and dispose of them in accordance with accepted standards. |
|  | Device classification type BF. |
| IP22 | Ingress protection rating. Leading digit (protection against ingress of solid foreign objects): 2 – protection against ingress of solid objects more than 12 mm in size; fingers or other objects with a maximum length of 80 mm, or solid objects. Second digit (protection against ingress of foreign liquids): 2 – protection against vertically falling water drops and objects when enclosure tilted up to 15° (normally positioned equipment). |
|  | Warning/Caution. |
|  | Consult user manual before use. |
|  | Authorized Representative in the European Union. |
|  | Manufacturer. |

3 SCOPE OF APPLICATION

Infrared thermometer AGU IHE3 is intended for measuring body temperature in temple area and in external ear canal in the home environment. Do not use the device for purposes other than intended. This thermometer can be used for children and adults.

4 COMPLETENESS

- 1 Infrared thermometer **AGU IHE3** – 1 pc.
- 2 Battery CR 2032 – 1 pc.
- 3 Storage case – 1 pc.
- 4 User manual – 1 pc.


5 BASIC SPECIFICATIONS

Specifications


| | |
|---------------------------|--|
| Type | Infrared thermometer |
| Model | AGU IHE3 |
| Power supply | 3 V CR 2032 |
| Temperature range | Measuring range: human body: 34°C~43°C (93.2°F~109.4°F) |
| Measuring accuracy | Laboratory accuracy: ±0.2°C (34°C~43°C)/±0.4°F (93.2°F~109.4°F) Calibration accuracy: ±0.2°C (35°C~42°C)/±0.4°F (95°F~107.6°F); out of range: ±0.3°C/±0.5°F Temperature measurement in external ear canal: ±0.2°C (34.0°C~43.0°C)/±0.4°F (93.2°F~109.4°F) Temperature measurement in the temple area: ±0.2°C (34.0°C~43.0°C)/±0.4°F (93.2°F~109.4°F) |
| Dimensions (LxWxH), mm | 125x35x34 |
| Weight, g | ~53 |
| Enclosure material | ABS plastic |
| Measuring time | About 1 second |
| Display | Liquid-crystal display, 4 digits with special icons and color display |
| Auto power-off | In 1 minute |
| Transport conditions | From -10°C to +40°C with relative humidity up to 75% |
| Storage conditions | From +10°C to +30°C with relative humidity up to 65% |

| | |
|----------------------|---|
| Operating conditions | 15°C~35°C (59°F~95°F) with relative humidity up to 15-95% |
| Memory | 10 previous temperature measurements |

6 PRODUCT DESCRIPTION

CONTROL PANEL AND DISPLAY

7 SAFE OPERATION INSTRUCTIONS

During early stages of illness with a rapid rise in temperature to high values, the effect of «white hyperthermia» can be observed – a condition which can be characterized by the vessel constriction, when the skin becomes pale and remains cold. In such cases, temperature measurements in temple area should not be carried out, because the temperature of the skin will be low. IN THIS CASE, YOU NEED TO MEASURE TEMPERATURE IN EXTERNAL EAR CANAL.

! ATTENTION

- The environmental temperature in the apartment may differ in various rooms. Before measuring temperature, the user and the thermometer should stay in the same room (premises) with constant environmental conditions for at least 30 minutes.
- Before starting measurement, please remove any dirt, hair or sweat from the measuring area.
- Do not measure temperature during or straight after nursing.
- The user should not take food, drink or perform physical exercises immediately prior to or during temperature measurement.
- Do not pocket the thermometer before use as it can be heated.
- Every time before starting measurement, please make sure that the measuring sensor lens is not contaminated. In case of contamination,

clean the measuring sensor lens, wait at least several minutes before further measurements in order to ensure the initial operating temperature.

When using this product, please observe the instructions listed below. Failure to follow these instructions may cause injury or affect the measurement accuracy:

- do not repair and remodel the thermometer;
- avoid contamination of the infrared sensor lens as it may result in incorrect readings;
- avoid direct finger contact with the lens;
- do not expose the thermometer to extreme temperature, very high humidity, or direct sunlight;
- avoid extreme shock or dropping the device;
- avoid measuring temperature within 30 minutes after exercise, bathing, or returning from outdoor;
- dispose of empty batteries at appropriate collection sites according to national or local regulations;
- it is ill-advised to disassemble the thermometer;
- use the thermometer for the intended purpose only.

8 PREPARATION FOR WORK

Remove the device from the package. Check the battery charge level.

When the battery gets weak, the battery warning symbol appears . It is still possible to measure temperature.


When the battery symbol flashes and «Lo» appears in the display, the battery must be replaced. If the battery is too flat, then the thermometer will switch off automatically.

NOTE

- It is recommended to remove the batteries if the thermometer will not be used for an extended period of time.
- Keep the battery away from children, do not expose to heat.
- Batteries must be disposed of in accordance with local environmental and institutional policies.
- Do not use rechargeable batteries.

BATTERY REPLACEMENT

- 1 Open the battery case by sliding off the cover on the back of the device.
- 2 Remove the used battery.
- 3 Replace with a lithium 3 V CR2032 battery. The larger part of the battery should be facing up. You should still be able to see the «+» sign when battery is installed.
- 4 Slide the cover back on and tighten the screw with a screwdriver. The device is ready for use.

Dispose of used batteries in accordance with the applicable legal regulations. Never dispose of batteries in the normal household waste.

9 OPERATION PRINCIPLE AND PROCEDURE

Bear in mind to keep the thermometer with the patient in the same room where the measurement is to be taken for at least 30 minutes before use.

MEASURING BODY TEMPERATURE IN THE EAR

! ATTENTION


- Some people produce different readings in their left and right ear. In order to record temperature changes, always measure a person's temperature in the same ear.
- The measurement must not be taken in an ear affected by inflammatory diseases, after possible ear injuries or during the warm-up period. In all of these cases, please consult your doctor.
- If you have been lying on one ear for some time, the temperature is slightly raised. Wait a little while or measure in the other ear.

- Accumulation of ear wax on the sensor can affect the measuring accuracy and result in the spread of infection among the thermometer users. It is extremely important for each measurement to observe the rules of hygiene and keep the sensor clean. To clean the sensor, follow the instructions in the section 10 «Cleaning». Every time before starting measurement, please make sure that the measuring sensor lens is not contaminated.
- After cleaning the measuring sensor lens with a disinfectant solution, wait at least several minutes before further measurements in order to ensure the initial operating temperature.


The thermometer takes temperature by measuring infrared heat generated by ear drum and surrounding tissues. Temperature measurement should be carried out by placing the tip in the external ear canal above the ear drum, and not over the surface of nearby anatomical areas. Failure to do so may result in excessive readings.

Body temperature measurement in the external ear canal ensures accurate readings due to the same temperature of blood flowing towards the ear drum and the brain's temperature regulation center – hypothalamus. Therefore, changes in body temperature are faster and more accurately recorded in the external ear canal than in any other body area.

- 1 Remove the front cover.


- 2 Press the Power/Measure button to turn on the device. Following a successful self-test, the device emits beep sound.
- 3 At the moment of switching on, all the icons are displayed on the display, after which the thermometer switches to the previous measurement mode (digital and light indication) and the device is ready for operation.

- 4 As the ear canal is slightly curved, you have to pull the ear slightly up and backwards before inserting the sensor tip. This is important so that the sensor tip can be pointed directly at the ear drum.
- 5 Insert the probe tip into ear canal carefully. Press the Power/Measure button.

- 6 The end of the measuring time is signaled with beep sound and the measured value appears on the display.
- 7 Read the recorded temperature from the LCD display. The backlight shows red for fever or green for normal.

! NOTE

≥37.5 °C: red backlight.

If the temperature measurement result is **37.5 °C or above**, the red backlight will light for 5 seconds.


<37.5°C: green backlight.

If the temperature measurement is **below 37.5 °C**, the green backlight will light for 5 seconds.



- 8 The thermometer will automatically turn off in 1 minute after you stop using it.

MEASURING BODY TEMPERATURE IN TEMPLE AREA

- 1 The cap should be in position on the probe nozzle.

- 2 Press the Power/Measure button. Following a successful self-test, the device emits beep sound.
- 3 At the moment of switching on, all the icons are displayed on the display, after which the thermometer switches to the previous measurement mode (digital and light indication) and the device is ready for operation.

- 4 Touch the probe gently on the skin of upper right of right eyebrow or the upper left of left eyebrow. If the skin area is covered with hair, sweat or dirt, please remove the obstacle beforehand to improve the reading accuracy.

- 5 Press the Power/Measure button and hold the probe flat on the area described above. You'll hear a beep and your temperature will come up on the screen.

- 6 Read the recorded temperature from the LCD display. The backlight shows red for fever or green for normal.

! NOTE

≥37.5 °C: red backlight.


If the temperature measurement result is **37.5 °C or above**, the red backlight will light for 5 seconds.

<37.5 °C: green backlight.

If the temperature measurement is **below 37.5 °C**, the green backlight will light for 5 seconds.

- 7 The thermometer will automatically turn off in 1 minute after you stop using it.

SWITCH BETWEEN CELSIUS (°C) OR DEGREES FAHRENHEIT (°F)


This thermometer can display result in either degrees Celsius (°C) or degrees Fahrenheit (°F).

When the device is on press the Power button and Memory button at the same time to toggle between °C and °F.

SETTING THE TIME

At the very first time of using this thermometer or when re-installing the battery, you may first set the clock time. Then, the ambient temperature measurement function will be automatically activated.

- 1 Open the battery cover, turn on the device, press the Time setting button.


- 2 When the hour numbers are blinking, click Power/Measurement button to adjust the correct hour.

- 3 Press the Time Setting button again to enter minute setting.

- 4 When the minute numbers are blinking, click Power/Measurement button to adjust the correct minute.



- 5 Press the Time setting button again to end setting.

When the Time function is activated, the icon  will be displayed.

CLOCK TIME AND AMBIENT TEMPERATURE

There are two display modes when activate clock function: **Clock time/Ambient temperature mode.**

After auto power off, the screen displays the time and temperature in an alternative way.

Note: using the device as a room thermometer and clock constitutes a continual measuring process, which significantly reduces the battery service life.

If user doesn't set the time in the first using or re-installing battery, there won't be display or time and ambient temperature.

MEMORY FUNCTION

You can recall up to 10 stored measurements in memory to share with your physician or trained healthcare professional.

- 1 Press the Memory button. The first reading displayed is the latest measurement stored in memory.
- 2 Continue to press the Memory button to view the next previously stored measurement.
- 3 Any new measurement will be recorded and the oldest memory deleted without you having to do anything.

When the Memory function is activated, the icon  will be displayed.

10 CLEANING



ATTENTION

Do not use chemically active agents to clean the device. Disinfecting agents can be used in the home.

Lens/Measurement sensor: in case of contamination of the infrared thermometer sensor lens, the readings begin to change immediately as something interferes with the measurement (usually the readings are too low). To eliminate this problem, it is enough to wipe the sensor with a soft lint-free cloth (or cotton swab) moistened with disinfectant solution. Then it is necessary to wipe the lens with a dry soft lint-free cloth. After the cleaning procedure,

the lens surface should be crystal clear. As the lens surface gets cool when the disinfectant solution evaporates, the thermometer should be left for several minutes at room temperature before use after cleaning.

Thermometer: clean with a soft, dry cloth. Do not use water to rinse the device.

11 TROUBLESHOOTING

If a problem occurs during operation, first of all, refer to the list of possible failures given below.

| Failure | Failure description | Remedial action |
|----------------|---|--|
| H _i | The temperature measured is higher than 43 °C (109.4 °F) | Operate the thermometer only between the specified temperature ranges. If necessary, clean the sensor tip. In the event of a repeated error message, contact your retailer or Customer Service |
| Lo | The temperature measured is lower than 34 °C (93.2 °F) | |
| Err | The operating temperature is not in the range 15~35 °C (59~95 °F) | Operate the device only at an ambient temperature within the specified values |

12 STORAGE, TRANSPORTATION AND OPERATION REQUIREMENTS

- The device shall be transported within the temperature range from -10 °C to +40 °C with relative humidity up to 75 %.
- The device can be operated within the temperature range from +15 °C to +35 °C with relative humidity up to 15–95 %.
- The device shall be stored within the temperature range from +10 °C to +30 °C with relative humidity up to 65 %.
- Do not expose to thermal shock.

ATTENTION

After transporting or storing at low temperatures, it is necessary to keep the device at room temperature for at least 2 hours before switching on.

13 DISPOSAL

The symbol on the product or its package means that this product does not fall under the category of domestic waste.

- Proper disposal of the device will prevent adverse environmental and human health effects.
- In order to protect the environment, the device must not be disposed of together with domestic (household) waste. Disposal shall be provided in accordance with local regulations.
- The device must be disposed of in accordance with the EU Directive 2012/19/EU WEEE (Waste Electrical and Electronic Equipment).

If you have any questions, please contact the local public utility responsible for waste disposal.

14 CERTIFICATION

This product conforms to the provisions of the EU Medical Device Directive (93/42/EEC). The following standards are applied to design and/or manufacture of the products:

- ISO 80601-2-56
Medical electrical equipment – Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement;
- IEC/EN 60601-1
Medical electrical equipment – Part 1: General requirement for safety;
- IEC/EN 60601-1-2
Medical electrical equipment – Part 2: Collateral standard: Electromagnetic compatibility – Requirements and tests.

This device is an adjusted mode clinical thermometer. The validated data for clinical accuracy in each adjusted mode are:

This device is an adjusted mode clinical thermometer. The validated data for clinical accuracy in each adjusted mode are:

- group A1:** $D_{cb} = -0,01 \text{ }^{\circ}\text{C}$, $L_A = 0,18$, $\sigma_r = \pm 0,08 \text{ }^{\circ}\text{C}$;
group A2: $D_{cb} = 0,06 \text{ }^{\circ}\text{C}$, $L_A = 0,22$, $\sigma_r = \pm 0,08 \text{ }^{\circ}\text{C}$;
group B: $D_{cb} = -0,01 \text{ }^{\circ}\text{C}$, $L_A = 0,20$, $\sigma_r = \pm 0,07 \text{ }^{\circ}\text{C}$;
group C: $D_{cb} = -0,01 \text{ }^{\circ}\text{C}$, $L_A = 0,18$, $\sigma_r = \pm 0,07 \text{ }^{\circ}\text{C}$.

D_{cb} – clinical bias, L_A – limits of agreement, σ_r – clinical repeatability.

ELECTROMAGNETIC COMPATIBILITY TABLES

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information.

- 1 Electromagnetic/Radio Frequency Interference (EMC/RFI): readings may be affected if the unit is operated within radio frequency electromagnetic field strength of approximately 3 volts per meter, but the performance of the instrument will not be permanently affected. Care should be taken to keep the Infrared thermometer **AGU IHE3** at least 6 inches away from R/C transmitters to avoid radio frequency interference.
- 2 Avoid keeping the Infrared thermometer **AGU IHE3** too close to objects that continuously generate high heat (like a hot plate) for long periods of time, which can cause overheating of the thermometer.

Table 1. For all ME EQUIPMENT and ME SYSTEMS

| Guidance and manufacturer's declaration – electromagnetic emissions | | |
|--|-------------------|---|
| Infrared thermometer AGU IHE3 is intended for use in the electromagnetic environment specified below. The customer or the user of the Infrared thermometer AGU IHE3 should assure that it is used in such environment | | |
| Emissions test | Compliance | Electromagnetic environment guidance |
| RF emissions CISPR 11 | Group 1 | Infrared thermometer AGU IHE3 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment |
| RF emissions CISPR 11 | Class B | Infrared thermometer AGU IHE3 is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes |
| Harmonic emissions IEC 61000-3-2 | Not applicable | |
| Voltage fluctuations/ Flicker emissions IEC 61000-3-3 | Not applicable | |

Table 2. For all ME EQUIPMENT and ME SYSTEMS

| Guidance and manufacturer's declaration – electromagnetic immunity | | | |
|---|---|--------------------------------|--|
| Infrared thermometer AGU IHE3 is intended for use in the electromagnetic environment specified below. The customer or the user of Infrared thermometer AGU IHE3 should assure that it is used in such environment | | | |
| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment guidance |
| Electrostatic discharge (ESD) IEC 61000-4-2 | ±6 kV (contact) ±8 kV (air) | ±6 kV (contact) ±8 kV (air) | Floor should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 % |
| Electrical fast transient/burst IEC 61000-4-4 | ±2 kV for power supply lines ±1 kV air for input/output lines | Not applicable | Mains power quality should be that of a typical commercial or hospital environment |
| Surge IEC 61000-4-5 | ±1 kV (differential mode) ±2 kV (common mode) | Not applicable | Mains power quality should be that of a typical commercial or hospital environment |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | <5 % U_T (>95 % dip in U_T for 0.5 cycle); 40 % U_T (60 % dip in U_T for 5 cycles); 70 % U_T (30 % dip in U_T for 25 cycles); 5 % U_T (>95 % dip in U_T for 5 sec) | Not applicable | Mains power quality should be that of a typical commercial or hospital environment. If the use of Infrared thermometer AGU IHE3 requires continued operation during power mains interruptions, it is recommended that Infrared thermometer AGU IHE3 be powered from an uninterruptible power supply or a battery |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 3 A/m | 3 A/m | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment |

NOTE: U_T – is the AC mains voltage prior to application of the test level.

Table 3. For EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

| Guidance and manufacturer's declaration – electromagnetic immunity | | | |
|--|--|----------------------------|---|
| <p>Infrared thermometer AGU IHE3 is intended for use in the electromagnetic environment specified below. The customer or the user of Infrared thermometer AGU IHE3 should assure that it is used in such environment</p> | | | |
| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment guidance |
| <p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p> | <p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p> | <p>3 Vrms</p> <p>3 V/m</p> | <p>Portable and mobile RF communications equipment should be used no closer to any part of AGU IHE3 infrared thermometer, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance: $d=1.2\sqrt{P}$, $d=1.2\sqrt{P}$ 80 MHz to 800 MHz, $d=2.3\sqrt{P}$ 800 MHz to 2.5 GHz,</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol: </p> |

NOTE 1: at 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: these guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a Field strengths from fixed transmitters, such as base stations from radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be adequately predicted by theoretical calculations. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which **Infrared thermometer AGU IHE3** is used exceeds the applicable RF compliance level above, **Infrared thermometer AGU IHE3** should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating **Infrared thermometer AGU IHE3**.
- b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.

Table 4. For EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

| Recommended separation distances between portable and mobile RF communications equipment and Infrared thermometer AGU IHE3 | | | |
|---|---|-------------------------------------|--------------------------------------|
| Infrared thermometer AGU IHE3 is intended for use in an electromagnetic environment where radiated RF distances are controlled. The customer or the user of Infrared thermometer AGU IHE3 can help to prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and Infrared thermometer AGU IHE3 as recommended below according to the maximum output power of the communications equipment | | | |
| Rated maximum output power of transmitter, W | Separation distance according to frequency of transmitter, m | | |
| | 150 kHz – 80 MHz d=1.2√P | 80 kHz – 800 MHz d=1.2√P | 800 kHz – 2.5 GHz d=2.3√P |
| 0.01 | 0.12 | 0.12 | 0.23 |
| 0.1 | 0.38 | 0.38 | 0.73 |
| 1 | 1.2 | 1.2 | 2.3 |
| 10 | 3.8 | 3.8 | 7.3 |
| 10 | 12 | 12 | 23 |

15 MANUFACTURER'S WARRANTY

The warranty period of this product is 24 months from the date of sale, providing that all operation conditions specified in this instruction manual are strictly observed. The warranty is valid subject to availability of a warranty card filled out by a duly authorized representative confirming the date of sale and also of the relevant receipt. The warranty does not cover the consumable parts, components and constituent elements, namely: battery, gift box, etc.

The warranty and free maintenance service are not provided in case:

- the device is used with violation of the requirements set by the operation manual;
- damage is caused by consumer's willful or wrong actions as a result of improper or negligent use;
- traces of mechanical effect (dents, cracks, chips, etc.), tampering, dismantling, traces of repairs outside the authorized technical service center, traces of ingress of moisture into the case or effect of corrosives, or of any other outside interference with the device design are observed on the device case, and also in other cases of violation by the consumer of the rules of storage, cleaning, transportation and maintenance of the device stipulated by the operation manual;
- penetration of oil, dust, insects, liquids (not intended for the use with devices) and other foreign matters inside the device.

Precisely follow the instructions to ensure reliable and long-term operation of the device.

In case of abnormal operation of the device, please contact your retailer or Customer Service.

The manufacturer reserves the right to make alterations in the design of the device.

ATTENTION

Do not perform unauthorized repair of the device – this will void the warranty.

The device alterations without the manufacturer's permission will void the warranty.

For repair and maintenance, please contact a specialized after-sales service (see agu-baby.com).

Date of manufacture  and importer information are specified on the individual package.

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AGU[®] is the registered trademark by Montex Swiss AG,
Medical Device Safety Service GmbH (MDSS),