



36.8 °C
M

37.5 °C
M



swissdesign

Giraffe

THE AGU NON-CONTACT THERMOMETER



non-contact



measurement in 1 second



high measuring precision



color and sound indicator



It is very important to get an accurate information on a child's temperature at first suspicion that he/she may be has a fever. But children do not like to stay still, and thus, they do not like getting their temperature taken. And during sleep, resting babies shouldn't be disturbed.

The AGU Non-contact Thermometer measures their temperature without touching the child. It takes only 3 to 5 seconds to do the reading with an accuracy of $\pm 0,2$ °C. It signals when the measurement taking has finished and it rings an alarm in case of fever. **Thermometer Giraffe AGU** can also measure surface and liquid temperatures.

agu-baby.com

MEASURING WITHOUT TOUCHING.